MARSDEN

USER MANUAL

BAS-200

Please take the time to read these instructions before starting to use the scales



Technical Specifications

- Working temperature: 10°C to 40°C
- Environmental humidity: <80%
- Max capacity: 20kg
- Graduation: 5g
- Power supply: 4 x AA batteries
- Background colour: Blue backlight
- Low power: Lo
- Error of indication: Err
- Unit change: kg or lb
- Memory function: Optionally save and recall weighing value
- Turn off automatically
- Overload: o_Ld

Keys' Function

ON/TARE	Press the key to power on the scale. To zero the scale, press the key again to clear the display and "zero" icon will show in the left corner of the LCD. To turn off the scales press and hold the key for about 3 seconds.
MODE	Weighing units exchange key. Press the key to select kg and lb in turn. An icon denoting kg or lb is on the right side of the LCD display.
MEM	When the display "0.000", press the key to recall the previous weight saved in the memory. A weight is saved by pressing the key when it is stable.

Operation

- Press "**ON/TARE**" key to start the scale until "**0.000**" appears.
- Place baby on the tray, the weight will flash after a few seconds. If the baby remains on the tray after the weight is shown, the weight will continue to flash for 5 minutes and then turn off automatically.
- If the baby is removed from the tray immediately after the weight has flashed, the scale will return to "**0.000**" and then turn off automatically within 5 minutes if not used again.
- To recall the saved weight, press the "**MEM**" key while the scale is in the "**0.000**" position.

Warnings

- The scales should be used on a hard flat surface.
- Replace batteries immediately when the scale displays "L o"
- Keep the scale clean
- Remove the batteries if the scale is not going to be used for a long period.
- Weigh the baby in the centre of the tray.
- Never overload the scale (23kg)



Head Office:	Unit 1
Sales, Service, Marketing &	Genesis Business Park
Accounts	Sh effi eld Road
	Rotherham
	S60 1DX
Telephone:	0845 13 0 7330
Fax:	0845 13 0 7440
Despatch, Production & Technical	lls it 7
Despacen, Froduction a recimicat	Unit 7
•	Centurion Business Park
•	
•	Centurion Business Park
•	Centurion Business Park Coggin Mill Way
Support	Centurion Business Park Coggin Mill Way Rotherham
Support Telephone:	Centurion Business Park Coggin Mill Way Rotherham S60 1FB