

# Marsden M-160 User Manual



## **Contents**

ntroduction	2
Product Specification	2
Safety Instructions	3
Explanation of Graphic Symbols	3
Assembling The Column	5
Before Use	5
Attaching the Height Measure	6
Power Supply and Low Battery	7
Operation: Basic Functions	8
Operation: Advanced Functions	10
Jsing the Scale with a Printer	13
Jsing the Scale with Bluetooth	14
EMC Guidance and Manufacturer's Declaration	15
Recommended Separation Distance	16
Error Messages	17
Manufacturer's Declaration of Conformity	17

## Introduction

Thank you for purchasing a professional gym weighing scale. Ideal for areas of heavy traffic and guaranteed accurate weight readings. Continue reading for instructions on how to use the device and company contact information.

# **Product Specification**

Model	M-160	
Capacity/Division	250kg x 100g	
Weight of Scale	Approximately 9kg	
Units of Measure	Kg/St	
Function Keys	ON/OFF, HOLD, TARE, BMI, UNIT, 0-9	
Stabilization Time	1-2 Seconds	
Operating Temperature	0 to 40°C	
	Rechargeable battery pack	
Power Supply	6 x AA batteries*	
	12V 1A AC Adaptor	
Indicator Display	2.5cm LCD display with 5 active digits	

#### **Safety Instructions**

Before putting the device into use, please read with care the information given in this user manual, which contains important instructions for proper installation, use and maintenance of the device.

Marsden and/or the manufacturer shall not be liable for damages arising from failure to heed the following instructions:

- When using electrical components under increased safety requirements, always comply with appropriate regulations.
- Inappropriate installation/use will render the warranty null and void.
- Ensure the voltage marked on the power supply unit matches your mains supply.
- This device is designed for use indoors only.
- Observe the permissible ambient temperatures for use.
- The device meets the requirements for electromagnetic compatibility. Do not exceed the maximum values specified in the applicable standards.
- Batteries should be kept away from small children. If swallowed, promptly seek urgent medical assistance.

If you have any problems with this scale, please contact Marsden/your local dealer/your service partner.

If a serious incident occurs in relation to this device, it should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

#### Cleaning

- We recommend using alcohol-based wipes or similar when cleaning the scale.
- Please do not use corrosive liquids, large amounts of water or high-pressure washers.
- Always disconnect the scale from the mains power supply before cleaning.

#### **Maintenance**

- The scale does not require any routine maintenance. However, we recommend checking the scale's accuracy at regular intervals. If any inaccuracies occur, please contact your local dealer or service partner.
- If you are in the UK, service contracts are available from Marsden to keep your scale accurate and reliable for longer.
   Call 01709 364296 for more information.

## Disposing of the Scale

- This product should not be treated as regular household waste but should be handed in to an electrical/electronic equipment recycling centre.
- You can obtain further details from your local council, your municipal waste disposal company or from where you
  purchased the product.
- Alternatively, you can return this product to Marsden we will recycle this free of charge,

#### Intended Use


- This scale is intended for use to determine the weight of patients, supported by professional personnel and in rooms intended for carrying out healthcare. The weighing value can be read after a stable weighing value has been obtained. Before use, the scale must be checked by an authorised person to ensure it's in a suitable condition.
- Device is intended to measure one subject at a time.

# **Explanation of Graphic Symbols**

$\triangle$	Caution, consult accompanying documents before use		Separate collection for waste of electrical and electronic equipment, in accordance with Directive 2002/96/EC
***	Manufacturer of medical device		Manufacturing year of medical device
	Carefully read user manual before installation and usage, and follow instructions for use.	*	Medical electrical equipment with Type B applied part
REF	Device catalogue number	EC REP	Authorized representative in the European Community
LOT	Manufacturer's batch or lot number	MD	Device is a medical device
SN	Serial number	UDI	Unique Device Identifier
<b>( (</b> 2460		Device conforms to 93/42/EEC as amended by Device Directive. Four digit number refers to N	
	D	Device complies with International Organizatio (Class III) requirements (verified models only)	n of Legal Metrology
<b>C€</b> M 190122		Device complies with EC directives (verified metals of the compliance with Directive automatic weighing instruments  19: Year in which conformity verification was pelabel was applied. (ex: 19=2019)  0122: Refers to Notified Body for metrology	re 2014/31/EU for non-
변 M21 0120		Device complies with UK Regulation.  M: Non-Automatic Weighing Instruments Regulation and the second	erformed and the CE

## Assembling the Column (M-160 only)

Assembling the M-160's column should ideally be carried out by two people.




Holding both the column and the base, ensuring all wiring is tucked inside the column housing, turn the scale over so that the column can be screwed to the base from the underside of the base.

Use the four screws provided, screwing them carefully into the four holes as shown in the drawing.

#### **Before Use**

- 1) Place the scale on a firm and level surface
- 2) Unwind the levelling feet until they make contact with the ground.



3) Ensure that all four of the levelling feet and the wheels are firmly on the ground and that the spirit level bubble is located in the center as shown below:




## Attaching the Height Measure (M-160 only)

The M-160 can be purchased with a height measure. If you did not purchase your M-160 with a height measure, it is also available separately (model HM-201M).



Using the four M5\*0.8\*9 screws provided, mount the two brackets onto the column of the M-160.



Using the four M5\*0.8\*25 screws provided, screw the black blocks onto the mounting brackets.



Using the two M6\*1.0\*10 screws provided, screw the HM- 201M to the black blocks.

To measure a person's height, retract the height rod and position the head stop accordingly.

For M-160 spare parts, call Marsden on 01709 364296.

## **Power Supply and Low Battery**

The indicator uses a rechargeable battery pack, a non-rechargeable battery pack, or can be powered from the mains via the AC adaptor.


Make sure the battery pack is installed in the battery box of the indicator. Alternatively, plug the AC adaptor (12V 2A) into the port on the side of the indicator.

## **Installing and Replacing the Battery Pack**

- 1. Take out the battery housing.
- 2. The rechargeable battery pack will slide into, or out of, the housing.




3. Check that the housing pin is connecting to the right point inside the indicator.




4. Place the housing back in the back of the indicator and close the battery housing cover.


## Switching on the Scale




## Switching off the Scale




## Setting the Scale to Zero




#### **Hold Function**




## **Setting Auto Hold Function (Optional)**



#### **Body Mass Index (BMI) Function**




## **Tare and Pre-Set Tare Functions**



## **Setting the Date**

Press HOLD for three seconds to access the time setting mode. The time period digit that is flashing can be changed by using the numeric keys. The time period to be edited is selected by pressing HOLD.

E.g., To input 25 December 2008, 8:00 a.m.:



## Using the Scale with a Printer

An optional Marsden external printer (model TP-2100) is available for all models. With the printer fitted, the patient's weight, height, and BMI result can be printed.

Once the person has been weighed and their BMI calculated, simply press PRINT to produce the following ticket:

GROSS WEIGHT 60.00kg

TARE WEIGHT 30.00kg

NET WEIGHT 30.00kg

PATIENT HEIGHT 100.0cm

PATIENT B.M.I 37. 5

29/12/2008 17:00


## **Connecting the TP-2100 Thermal Printer**



Plug the cable to the printer, and then connect its 9D connector to the indicator.

## Using the Scale with Wi-Fi/Bluetooth

If your scale has device connectivity, the universal wireless symbol will be on the main indicator display. If your scale does have Wi-Fi or Bluetooth connectivity, we do not supply the software to capture the data from the weighing scale. We do however provide the protocols for you to implement the devices into your own software.



## **Bluetooth Connection**

ROFF	Long press ZERO for three seconds to enter the Setting mode and then display the A-OFF menu.
P InEF	Press TARE twice, and then press HOLD once to enter the Bluetooth setting mode.
On ← → OFF	Using the HOLD button, select "ON" (enable) or "OFF" (disable). Press TARE to confirm the setting. Note: Disabling the Bluetooth function when not in use will reduce battery power consumption.
P InEF	Display the "bluEt" menu. Press TARE once.
	Press HOLD to return to normal mode.
End	Search for the scale in your computer or device's Bluetooth settings (procedure may vary depending on device or system)
LIU	The scale will appear on the Bluetooth device list as "MARSDEN BT".
	Connect your device to "MARSDEN BT", and the scale is ready to transmit data wirelessly via Bluetooth.

## **Wi-Fi Connection**

H.F.	Turn the scale on and press the TARE key for 3 seconds to enter into settings. Press the HOLD key several times until the display shows the information on the left.  When the display shows the symbols shown on the left, press the TARE key to see the status (ON/OFF).  If the display shows OFF, press the HOLD key once and it will change to ON. Then press the TARE key to confirm the setting.
End	Press the HOLD key several times until END shows on the display. When END is on the display, press the TARE keypad to enter normal weighing mode.
Spella, kaj jeforje audu Zo	If the highlighted triangular sign is not visible, Wi-Fi is turned off. If the triangle is solid, the device is connecting. When the triangle is blinking, the device is successfully connected.

#### **EMC Guidance and Manufacturer's Declaration**

Guidance and manufacturer's declaration – electromagnet emissions.

The M-160 is intended for use in the electromagnetic environment specified below. The customer or user of this scale should ensure that it is used in such environment.

Emission Test	Compliance	Electromagnetic environment- guidance
RF emissions CISPR 11	Group 1	This scale uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This scale is suitable for use in all establishments, including domestic
Harmonic emissions IEC 61000-3-2	Class A	establishments and those directly connected to the public low-voltage power supply network that supplies
Voltage fluctuations/flicker emissions IEC 61000-3-3	Compliance	buildings used for domestic purposes.

## Guidance and manufacturer's declaration - electromagnetic immunity.

The M-160 is intended for use in the electromagnetic environment specified below. The customer or the user of this scale should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, cement or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines +1 kV for input/output lines	± 2 kV for power supply lines not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV line(s) to line(s) ± 2 kV line(s) to earth	1 kV differential mode not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5s	<5% UT (95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5s	Mains power quality should be that of a typical commercial or hospital environment. If the user of this scale requires continued operation during power mains interruptions, it is recommended that this scale is powered from an uninterruptable power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	The scale's power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note UT is the A.C mains voltage prior to application of the test level.			

## Guidance and manufacturer's declaration – electromagnetic immunity.

This scale is intended for use in the electromagnetic environment specified below. The customer or the user of the scale should ensure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment- guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 KHx to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the scale including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

Recommended separation distance:

 $d = 1,2 \sqrt{P}$ 

 $d = 1.2 \sqrt{P} 80MHz$  to 800 MHz

 $d = 2.3 \sqrt{P} 800MHz to 2.5 GHz$ 

Where *P* is the maximum output power rating of the transmitter in watts (w) according to the transmitter manufacturer and *d* is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.

Interference may occur in the vicinity of equipment marked with the following symbol:

(('<u>`</u>'))

Radiated RF IEC 61000-4-3 3 V/m 80 MHz to 2,5 GHz 3 V/m

NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- A) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the scale is used exceeds the application RF compliance level above, the scale should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the scale.
- B) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distance between portable and mobile RF communications equipment and the M-160.

This scale is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the scale can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the scale as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output	Separation distance according to frequency of transmitter m		
power of transmitter.	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz
W	$d = 1,2\sqrt{P}$	d = 1,2√ <i>P</i>	
		·	d = 2,3√ <i>P</i>
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output rating of the transmitter in watts (w) according to the transmitter manufacturer.

NOTE1) At 80 MHz and 800 MHz, the separation distance for the high frequency range applies.

NOTE2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## **Error Messages**

Low Battery The scale's alkaline AA type batteries are flat; please replace the batteries.	Lo
Overload This indicates that the scale's load sensor(s) have been overloaded. Reduce the loading and retry.	Err
The signal from the load cells is too high. Please remove any weight from the scale and try to power on again. If the scale continues to show the error message, it indicates a fault with the electronics or wiring.  The signal from the load cells is too low. Please remove any weight from the scale and try again. If	Err.H Err.L
the scale continues to show the error message, it indicates a fault with the electronics or wiring.  High/Low Zero Count  1. The scale is above its zero range. Please remove any weight from the scale and power on again. If the scale continues to show the error message, it	00000
<ul> <li>indicates a fault with the electronics.</li> <li>The scale is below its zero range. Check there is nothing jammed underneath the scale and power on again. If the scale continues to show the error message, it indicates a fault with the electronics.</li> </ul>	00000
EEPROM Error This indicates there is a fault with the scale's software and is normally caused by a fault with the load cell or wiring. Contact your local service representative.	Err.P

EU Authorized Representative:	C REP Obelis s.a.  Bd Général Wahis, 53 B-1030 Brussels Belgium
Distributor:	MARSDEN  Marsden Weighing Machine Group Ltd, Unit 1, Genesis Business Park, Sheffield Road, Rotherham, UK, S60 1DX
EU Importer:	MARSDEN  Marsden Weighing Machine Group Europe Ltd, The Black Church, St. Mary's Place, Dublin 7, Dublin, Ireland, D07 P4AX
Manufactured by:	Charder Electronic Co., Ltd. No.103, Guozhong Rd., Dali Dist., Taichung City 41262 ,Taiwan (R.O.C.)

# **MARSDEN**

Unit 1, Genesis Business Park, Sheffield Road, Rotherham, S60 1DX

Telephone: + 44 (0) 1709 364296

Version 2.0