# MARSDEN

# USER MANUAL

Please take time to read these instructions before starting to use the scale

# Contents

Introduction	3
Product Specification	3
Safety Instructions	4
Explanation of Graphic Symbols	5
Power Supply and Low Battery	6
Installing & Replacing the Battery Pack	6
Your M-900	7
Setting up the Scale	8
Operation: Basic Functions	10
Operation: Advanced Functions	12
Using the Scale with a Printer	16
Using the Scale with Bluetooth	17
EMC Guidance and Manufacturer's Declaration	18
Recommended Separation Distance	20
Error Messages	21
Manufacturer's Declaration of Conformity	22

Thank you for purchasing this Marsden professional medical scale. This is a precision Class III Weighing Instrument and considerate use will result in many years of accurate weighing.

The scale has a maximum load capacity of 600kg which must not be exceeded.

Model	M-900
Accuracy Class	Class III
Capacity/Division	600kg x 200g
Weight of each beam	Approximately 11.5kg
Units of Measure	Кд
Function Keys	ON/OFF, ZERO, HOLD, TARE, BMI, UNIT, 0-9
Stabilisation Time	1-2 Seconds
Operating Temperature	5 to 35°C
Power Supply	Rechargeable battery pack (7.2V 2000mA) 12V 1A switching adaptor (UE24WV- 120100SPA & UE24WB-120100SPA)
Indicator Display	2.5cm LCD display with 5 <sup>1</sup> / <sub>2</sub> active digits
Dimensions (weighing area)	1000mm x 160mm per beam

# **Product Specification**

# Safety Instructions

Before putting the device into use, please read with care the information given in this user manual, which contains important instructions for proper installation, use and maintenance of the device.

Marsden/the manufacturer shall not be liable for damages arising from failure to heed the following instructions:

- When using electrical components under increased safety requirements, always comply with appropriate regulations.
- Inappropriate installation/use will render the warranty null and void.
- Ensure the voltage marked on the power supply unit matches your mains supply.
- This device is designed for use indoors.
- Observe the permissible ambient temperatures for use.
- The device meets the requirements for electromagnetic capability. Do not exceed the maximum values specified in the applicable standards.
- Batteries should be kept away from small children. If swallowed, promptly seek medical assistance.

If you have any problems, contact Marsden/your local dealer/your service partner.

# Cleaning

- We recommend using alcohol-based wipes or similar when cleaning the scale.
- Please do not use corrosive liquids, large amounts of water or high pressure washers.
- Always disconnect the scale from the mains power supply before cleaning.

## Maintenance

- The scale does not require any routine maintenance. However, we recommend checking the scale's accuracy at regular intervals. If any inaccuracies occur, please contact your local dealer or service partner.
- To keep your scale calibrated, Marsden service contracts are available in the UK. Call 01079 364296 to find out more.

# Disposing of the Scale

- This product should not be treated as regular household waste, but should be handed in to an electrical/electronic equipment recycling centre.
- You can obtain further details from your local council, your municipal waste disposal company or from where you purchased the product.

# **Explanation of Graphic Symbols**

## SN-21300100





Designation of the serial number of every device. (Number as an example)

"Please note the accompanying documents" or "Observe operating instructions"

Identification of manufacturer of medical product including address.


Charder Electronic Co. Ltd No.103 Guozhong Rd, Dali Dist, Taichung City 412, Taiwan (R.O.C)

> Ť T

Type B applied part.

Dispose of old appliances separately from your household waste. This product must be disposed of at a communal collection point.

Carefully read this operation manual before setup and commissioning, even if you are already familiar with Marsden scales.



Transport and storage temperature limit indicating the upper and lower limit (transport and storage temperature on packaging).

# Power Supply & Low Battery


The indicator uses a rechargeable battery pack, or can be powered from the mains via the AC adaptor.

Make sure the battery pack is installed in the battery box of the indicator. Alternatively, plug the AC adaptor (12V 1A) into the port on the side of the indicator.




# Installing & Replacing the Battery Pack

- 1. Take out the battery housing.
- 2. The rechargeable battery pack will slide into, or out of, the housing.



3. Check that the housing pin is connecting to the right point inside the indicator.



4. Place the housing back in the back of the indicator, and close the battery housing cover.


## Your M-900




#	Part Name
1	Weigh Beam – 1
2	Weigh Beam – 2
3	Interconnecting Cable
4	Indicator
5	AC Adaptor

# Setting up the Scale


1. The bed scale has been labeled for easy assembly. Connect the weigh beam 1 & 2 using the interconnecting cable (3) and then connect the indicator (4) with weigh beam (2) as below:




2. Using the handle on each weigh beam...



3. ...place the weigh beams under the bed next to the castors.



- 4. Ensure that the beams are level by checking the bubble indicator. If the bubble is in the centre of the circle, the scale is levelled correctly. All four feet of each beam must touch the surface to ensure that M-900 does not move whilst in use.
- 5. Switch scale on, check the indicator display shows zero and push the bed forward and onto the weigh beams. We strongly recommend this is carried out by two people.






You must make sure that the connecting cable is not caught under the feet of the beams, otherwise faulty measurements may result.

# **Operation: Basic Functions**

## Switching on the Scale








Once the scale has tested all segments of the display, it will go into weighing mode and should show 0.0 on the display.

### Switching off the Scale




Press ON/OFF when the scale is turned on. The scale will now power down.

## Setting the Scale to Zero



If for any reason the scale shows a reading other than 0.0 it can be reset to zero.


Press ZERO once.




The scale will return to 0.0.

# **Operation: Advanced Functions**

## Hold Function










Press HOLD once.

Push the bed onto the beams.

After a few seconds the scale will lock on the weight. When the patient leaves the scale, or the bed is removed from the scale, the weight will remain on the display.













In normal mode, press BMI to enter BMI mode.

The display will show the last height entered and the extreme left digit will flash. Enter the height by using the numeric keys. Press ZERO to confirm the height. (NB: There will always be an active flashing digit in the height display, unless HOLD is pressed).

Weigh the patient as normal. The display will show the weight, height and BMI value. At this time, the weight and height can be freely changed, and the BMI value will be automatically calculated according to the changed weight and height.




Press BMI to return to normal weighing mode.




## 

Press TARE for three seconds to enter Preset Tare setting mode.

Press TARE again once. The display will show the last preset tare entered and the extreme left digit will flash. Enter the preset tare value by using the numeric keys, then press the TARE key again to confirm the value.



The display will show the figure you entered, deducted from 0.0kg. Press the ZERO key to return to normal weighing mode.



To use the Tare function, add the item you wish to tare off to the scale, and press the TARE key. The display will show zero, and then a minus number when the item is removed from the scale. Press HOLD for three seconds to access the time setting mode. The time period digit that is flashing can be changed by entering the appropriate number from the numeric key pad. The time period to be edited is selected by pressing HOLD.

E.g. To input 25 December 2008, 8:00 a.m.:





Enter the year. Press HOLD to confirm and access the date editing field.

Enter the date. E.g. "12.25" for December 25th. Press HOLD to confirm and access the time clock editing field.



Enter the time (24 hour clock only).

≥008	225 🔿	0800
------	-------	------

Press HOLD and the display shows: YYYY $\rightarrow$ MM.DD $\rightarrow$ HH:SS

0

Press HOLD to return to normal weighing mode.

# Using the Scale with a Printer

An optional Marsden external thermal printer (model TP-2100) is available for this scale. When the printer is fitted, the patient's weight, height, and BMI result can be printed.

Once the person has been weighed and their BMI calculated, simply press PRINT to produce the following ticket:

GROSS WEIGHT	60.00kg
TARE WEIGHT	30.00kg
NET WEIGHT	30.00kg
PATIENT HEIGHT	100.0cm
PATIENT B.M.I	37.5
29/12/2008	17:00


Connecting the TP-2100 Thermal Printer




Plug the cable to the printer, and then connect its 9D connector to the indicator.

# Using the Scale with Bluetooth

If your scale has Bluetooth connectivity, the universal wireless symbol will be on the main indicator display.



#### **Bluetooth Connection**



Long press the ZERO key for three seconds to enter the Setting mode and then display the A-OFF menu.

Press the TARE key twice, and then press HOLD once to enter the Bluetooth setting mode.

Using the HOLD key, select "ON" (enable) or "OFF" (disable). Press the TARE key to confirm the setting.

Note: Disabling the Bluetooth function when not in use will reduce battery power consumption.

Display the "bluEt" menu. Press the TARE key once.

Press the HOLD key to return to normal mode.

Search for the scale in your computer or device's Bluetooth settings (procedure may vary depending on device or system)

The scale will appear on the Bluetooth device list as **"MARSDEN BT".** 

Connect your device to "MARSDEN BT", and the scale is ready to transmit data wirelessly via Bluetooth.

# EMC Guidance and Manufacturer's Declaration

#### Guidance and manufacturer's declaration – electromagnet emissions.

The M-900 is intended for use in the electromagnetic environment specified below. The customer or user of this scale should ensure that it is used in such environment.

Emission Test	Compliance	Electromagnetic environment- guidance
RF emissions CISPR 11	Group 1	This scale uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This scale is suitable for use in all establishments, including domestic
Harmonic emissions IEC 61000-3-2	Class A	establishments and those directly connected to the public low-
Voltage fluctuations/flicker emissions IEC 61000-3-3	Compliance	voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration – electromagnetic immunity.

The M-900 is intended for use in the electromagnetic environment specified below. The customer or the user of this scale should ensure that it is used in such an environment.

IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
±6 kV contact ± 8 kV air	<ul> <li>± 6 kV contact</li> <li>± 8 kV air</li> </ul>	Floors should be wood, cement or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
± 2 kV for power supply lines +1 kV for input/output lines	± 2 kV for power supply lines not applicable	Mains power quality should be that of a typical commercial or hospital environment.
$\pm$ 1kV line(s) to line(s) $\pm$ 2 kV line(s) to earth	± 1 kV differential mode not applicable	Mains power quality should be that of a typical commercial or hospital environment.
<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5%	<5% UT (95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5%	Mains power quality should be that of a typical commercial or hospital environment. If the user of this scale requires continued operation during power mains interruptions, it is
	<pre>± 6 kV contact ± 8 kV air</pre> ± 2 kV for power supply lines +1 kV for input/output lines ± 1kV line(s) to line(s) ± 2 kV line(s) to earth <5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip	± 6 kV contact± 6 kV contact± 8 kV air± 6 kV contact± 8 kV air± 8 kV air± 2 kV for power supply lines +1 kV for input/output lines± 2 kV for power supply lines not applicable± 1 kV line(s) to line(s) ± 2 kV line(s) to earth± 1 kV differential mode not applicable<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5%

	5s	5s	recommended that this scale is powered from an uninterruptable power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	The scale's power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note UT is the A.C mains voltage prior to application of the test level.			

### Guidance and manufacturer's declaration – electromagnetic immunity.

This scale is intended for use in the electromagnetic environment specified below. The customer or the user of the scale should ensure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 KHx to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the scale including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
			Recommended separation distance: $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ 800 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz Where <i>P</i> is the maximum output power rating of the
			transmitter in watts (w) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by
			an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked

Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	with the following symbol:
NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies.			

NOTE: At so winz and boo winz, the higher frequency range applies. NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- A) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the scale is used exceeds the application RF compliance level above, the scale should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the scale.
- B) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

# Recommended separation distance between portable and mobile RF communications equipment and the M-900.

This scale is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the scale can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the scale as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output	Separation distance according to frequency of transmitter m		
power of transmitter	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz
W	d = 1,2√ <i>P</i>	d = 1,2√P	d = 2,3√ <i>P</i>
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output rating of the transmitter in watts (w) according to the transmitter manufacturer.

NOTE1) At 80 MHz and 800 MHz, the separation distance for the high frequency range applies. NOTE2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

# Error Messages

Low Battery The scale's alkaline AA type batteries are flat; please replace the batteries.	Lo
Overload This indicates that the scale's load sensor(s) have been overloaded. Reduce the loading and retry.	Err
<ul> <li>Counting Error</li> <li>1. The signal from the load cells is too high. Please remove any weight from the scale and try to power on again. If the scale continues to show the error message, it indicates a fault with the electronics or wiring.</li> <li>2. The signal from the load cells is too low. Please remove any weight from the scale and try again. If the scale continues to show the error message, it indicates a fault with the electronics or wiring.</li> </ul>	Err.H Err.L
High/Low Zero Count 1. The scale is above its zero range. Please remove any weight from the scale and power on again. If the scale continues to show the error message, it indicates a fault with the	0000
electronics. 2. The scale is below its zero range. Check there is nothing jammed underneath the scale and power on again. If the scale continues to show the error message, it indicates a fault with the electronics.	0000
EEPROM Error This indicates there is a fault with the scale's software and is normally caused by a fault with the load cell or wiring. Contact your local service representative.	Err.P

# Manufacturer's Declaration of Conformity



## Manufactured by:



Charder Electronic Co., Ltd. No.103, Guozhong Rd., Dali Dist., Taichung City 412, Taiwan (R.O.C.)



## **EU Declaration of Conformity**

The Non-Automatic Weighing Instrument

III

Manufacturer	Charder Electronic Co., Ltd
Model	M-900
EC Type Approval Certificate No.	T7616

The Metrological Aspects of Non-Automatic Weighing Instruments

EN45501:2015 (module D)	Notified Body Number – 0126
EN45501:1992 (module B)	Notified Body Number – 0122

The non-automatic weighing instrument corresponds to the production model described in the EC Type Approval Certificate and requirements of the following EC Directives:

2014/31/EU	Non-Automatic Weighing Instruments Directive
2014/30/EU	Electromagnetic Compatibility Directive
2014/35/EU	Low Voltage Directive

The applicable harmonized standards are:

EN45501:2015	The Metrological Aspects of Non-Automatic Weighing Machines
EN 61000-3-2:2014	Electromagnetic compatibility (EMC) - Part 3-2: Limits - Limits for harmonic
	current emissions (equipment input current ≤ 16 A per phase)
EN 61000-3-3:2013	Electromagnetic compatibility (EMC) - Part 3-3: Limits - Limitation of
	voltage changes, voltage fluctuations and flicker in public low-voltage
	supply systems, for equipment with rated current <= 16 A per phase and
	not subject to conditional connection
EN 60950-1:2006	Information technology equipment - Safety - Part 1: General requirements

This declaration of conformity is issued under the sole responsibility of the manufacturer.

Date: August 25, 2017

Signature:

gela In

Name: Angela Lu Position: Measuring Management Rep. Place: Taichung, Taiwan

Manufacturer: Charder Electronic Co., Ltd.

Address: NO.103, Guozhong Rd., Dali Dist., Taichung City 412, Taiwan (R.O.C.)



Accuracy Assured

Tel: 01709 364296 / 0800 169 2775 Fax: 01709 364293 E-mail: sales@marsdengroup.co.uk

#### Manufacturing and Distribution:

Unit 7, Centurion Business Park, Coggin Mill Way, Rotherham, S60 1FB Head Office: Unit 1, Genesis Business Park, Sheffield Road, Rotherham S60 1DX

IN-8587 80001B