

M-320 Ultra Lightweight Baby Scale

- Class III
- Accurate to 20kg x 10g / 50kg x 20g<20kg>50g
- Perfect for use in the community or the surgery
- Powered by 4 x AA batteries
- 100 hours of continuous use
- Auto power-off function helps to save battery power
- Ultra lightweight and highly portable
- Optional carry case available to aid portability
- Easy to read display and built-in spirit level
- Hold and Tare functions to enable accurate and easy weighing
- Made from Acrylonitrile Butadiene Styrene (ABS)
- Covered by the Marsden 8 year warranty



Versatile weighing

The M-320 is a versatile baby scale which can be used to weigh both babies and toddlers with optional capacity of either 20kg or 50kg. Slide off the tray to reveal a platform for toddlers to stand on. The scale is constructed from Acrylonitrile Butadiene Styrene to ensure that the structure is strong and the scale will be reliable and can be used for a long duration of time.

Ideal for the community

We understand that babies and toddlers are weighed at home, just as often as they are in a medical environment and because of this the M-320 is ultra lightweight and portable. This device does have an optional RK-400 Ruck sack, CC-400 Carry Case available to make transporting the scale even easier for community midwives and health visitors. It also adds a layer of protection when in transit.

Easy to use

With four simple buttons and a large easy to read display, the M-320 could not be simpler to use and this is why it is one of our best-selling products. The On/Off button ensures that you have full control over the device. Hold enables you to gain weights even for the most active babies. The tare keys allow you deduct the weight of the sling or any clothing. The recall memory button is helpful for regular client visits.

Capacity:	20kg / 50kg	Power Supply:	4 x AA Batteries
Graduations:	20kg x 10g / 50kg x 20g<20kg>50g	Battery Life:	100 hours of continuous use
Functions:	On/Off, Zero/Tare/Send, Recall/Memory, Hold/Milk Intake	Approvals:	Class III
Dimensions:	560mm x 308mm x 140mm	Optional Extras:	RK-400 Rucksack/CC-400 Carry Case
Weight:	20kg -.2.56kg / 50kg - 2.69kg	Warranty Period:	8 Years