

M100 Assess

Premium column scale supporting MUST, WHtR, with integrated connectivity for seamless synchronisation with EMR systems.

In partnership with The British Association for Parenteral and Enteral Nutrition.

- Class III Approved
- Accurate to 300kg x 100g
- Integrated automatic digital height measurement
- Built-in Wi-Fi and Bluetooth connectivity
- Large, stable platform for patient comfort and safety
- Optional printer for immediate documentation
- Touchscreen indicator with intuitive interface
- Patient ID via on-screen entry or barcode scanner



Streamlined Workflows

Can capture multiple data inputs during the same weighing event instead of weighing and then using separate charts, calculators or manual look-up tables, resulting in a MUST score, WHtR, and BMI calculations within one indicator.

Built-in Connectivity

Ready for connected healthcare, the built-in Wi-Fi and Bluetooth allow the scale to integrate easily into hospital and clinic networks, supporting digital records and modern workflows without additional add-ons.

Optional Barcode Scanner and Printer

Optional barcode scanner simplifies and speeds up data input for assigning a patient ID.

Optional printer allows calculations to be printed at the point of measurement for records or for the patient to take away.

Capacity:	300kg	Power Supply:	Rechargeable battery pack
Graduations:	300kg x 100g	Battery Life:	Up to 40 hours continuous use
Functions:	Tare/Hold/Auto Hold/Auto Power-off	Approvals:	Class III
Dimensions:	460mm x 440mm x 1290mm	Optional Extras:	Printer, Barcode Scanner
Weight:	16kg	Warranty Period:	8 years