

## M100 Sync

Column scale with integrated connectivity for seamless synchronisation with EMR systems.

- Class III Approved
- Accurate to 300kg x 100g
- Integrated automatic digital height measurement
- Built-in Wi-Fi and Bluetooth connectivity
- Large, stable platform for patient comfort and safety
- Optional printer for immediate documentation
- Touchscreen indicator with intuitive interface
- Patient ID via on-screen entry or barcode scanner



### Built-in Connectivity

Ready for connected healthcare, the built-in Wi-Fi and Bluetooth allow the scale to integrate easily into hospital and clinic networks, supporting digital records and modern workflows without additional add-ons.

### Automatic Height Measurement

The integrated height measure, with a maximum height 202cm, directly connects to the indicator, to provide fast, consistent results without additional equipment or height recording, saving time and reducing manual handling and human error.

### Large, Stable Platform

An oversized base plate provides confident footing for a wide range of patients, improving comfort and safety during measurement.

Capacity:	300kg	Power Supply:	Rechargeable battery pack
Graduations:	300kg x 100g	Battery Life:	Up to 40 hours continuous use
Functions:	Tare/Hold/Auto Hold/Auto Power-off	Approvals:	Class III
Dimensions:	460mm x 440mm x 1290mm	Optional Extras:	Printer, Barcode Scanner
Weight:	16kg	Warranty Period:	8 years