

BENCH SCALES

DS-673SS

DIGI TRADE APPROVED WATERPROOF BENCH SCALES

Features

- EC approved Class III bench scale
- Accuracy - 1/7500 approved
- IP68 certified
- Stainless steel construction
- Rechargeable battery
- Ideal for harsh environments
- AC adaptor
- Long battery life
- Large easy to read LED display
- Tare
- Zero



Product Code	Capacity	Class III Graduation	IP Rating	Mains Rechargeable	Pan Dimensions
DS-673SS-6	6kg	1g	68	Mains/Battery	250 x 290mm
DS-673SS-15	15kg	2g	68	Mains/Battery	250 x 290mm

