

# MARSDEN

## M-425

ITEM NO	DESCRIPTION	PART NO	PRICE £
1	OVERLAY	MP20000141	12.00
2	TILE	BLACK	12.00
3	DISPLAY BOARD	90016001471	25.00
4	LOAD CELL 1022N (150KG)	ME93000191	55.00
5	A/D BOARD	CH-8371	75.00
6	FEET	MP00600192	4.00
7	BATTERY BOX	MP00001691	4.00
8	BATTERY COVER	AE-0162	4.00
9	SLIDING SWITCH ON/OFF	AE-0992	4.00
10	AC JACK SOCKET	SCD-523	7.00
11	VR BOARD ASSEMBLY	90016001475	12.00
12	ADAPTOR 12V 300mA (TIP +)	ILD41B-1200300	36.00
13	REMOTE DISPLAY		40.00

Sales: 0845 130 7330  
[sales@marsdengroup.co.uk](mailto:sales@marsdengroup.co.uk)

Service: 01709 364 296  
[service@marsdengroup.co.uk](mailto:service@marsdengroup.co.uk)

Website:  
[marsden-weighing.co.uk](http://marsden-weighing.co.uk)