

MARSDEN

M-125

Marsden Column Scale



Class III Approved, 250kg capacity medical scale

100g accuracy - perfect for wards, clinics and outpatients

Simple indicator with large, clear 30mm display

Optional height measure

www.marsden-weighing.co.uk



Marsden M-125: Our entry level column scale, upgraded

Marsden's M-120 column scale has changed. Now available with our **brand new I-300 indicator**, this scale is still robust, reliable and long-lasting. But functions are **easier to access**, weight readings are **clearer and more accurate** - and **data can be transferred directly from the scale to a PC**.

Super accuracy for medical purposes

The M-120 was accurate to 200g. The M-125 is twice as accurate with **100g graduations**, yet still retaining the **250kg capacity** needed to allow weighing of **almost any patient**.

100g accuracy makes the M-125 suitable for **assisting in medical diagnosis**, as well as **monitoring weight for critical treatment**.

Easy access features to save time

The indicator for the scale has been **completely redesigned**, which has allowed us to make the **display larger** and features more **easily accessible**. Fewer buttons pressed means less time spent weighing.

Use the M-125 to measure weight, **BMI** and **BSA** - with the optional height measure - height as well. **Hold** stabilises the weight on the screen. **Tare** will remove unwanted weight, like clothing.

The Body Surface Area (BSA) calculation is, we believe, **a world first for weighing scales**. BSA is used to calculate drug dosages, primarily in hospitals, and on the M-125 it's as simple to calculate as BMI. The M-125 calculates BSA using the **Du Bois formula**.

www.marsden-weighing.co.uk





Solid, robust construction for years of reliable weighing

The M-125 has been engineered to **last year after year**. This is not a weighing scale you can expect to replace after three years. Used properly, an M-125 will outlive similar scales from almost any manufacturer.

This is because it has a **reinforced aluminium base** built for regular use, a rigid column and **wheels in the base of the column for easier transportation**. Plus, **6x AA batteries will last up to 180 hours before they need replacing** - that's around **9000 weigh-ins**.

Key features

- Class III Approved
- Mains operated or batteries
- Large, bright display
- Battery status indicator
- Wheels on base for transportation
- BMI/BSA calculation
- Hold/Tare
- Optional height measure



Accuracy assured for GPs, hospitals and health centres

The 100g accuracy of the M-125 means it's perfect for medical weighing for GPs, hospitals and health centres.

It could also be a weighing solution in waiting rooms - patients could use the M-125 with ease, since just five buttons make it **easy to use when you simply want to see your weight**.

Specification

- Class: III
- Capacity: 250kg
- Graduations: 100g
- Power: Mains adaptor or 6 x AA batteries
- Base dimensions: 310mm x 90mm x 310mm
- Weight of Scale: 9kg

Marsden M-125: Key features at-a-glance



Tel: 01709 364296 / 0800 169 2775

Fax: 01709 364293

E-mail: sales@marsdengroup.co.uk

www.marsden-weighing.co.uk

 MarsdenWeighing

 @MarsdenWeighing

 www.youtube.com/MarsdenWeighing

 <http://pinterest.com/marsdenweighing>

MARSDEN