

# MARSDEN

## V-20

Marsden Portable Vet Scale






20kg capacity animal scale with carry handles

Powered by a single 9v battery

Weighs in kg or lb

Suitable for weighing smaller pets

[www.marsden-weighing.co.uk](http://www.marsden-weighing.co.uk)



## Marsden V-20: A premium, portable vet scale ideal for surgery or home use

The V-20 is a **small and portable veterinary scale** which is **ideal for use in both the home for domestic pets and in veterinary surgeries** for a wider range of animals. Using the V-20 is great way to keep track of an animal's overall health.

### Weigh your pet - anywhere

The V-20 is a **20kg superbly designed vet scale**, built for weighing pets anywhere. It's **easy to use**, and is **powered by a 9v battery**.

This digital weighing scale is **lightweight at only 3.5kg** and there is an **optional carry case** which makes the scale **even easier to move and store**. Plus, the scale is **ultra slim** and **easy to clean**.

### Quick and easy to use

The V-20 has a **Tare function** so any unwanted weight can be deducted from the overall weight reading. An accurate weight reading is also always guaranteed thanks to the **Hold function**, which displays a stable weight even when the animal is moving on the scale.

The weight reading, **accurate to the nearest 10g**, is shown clearly on the display.



## Use absolutely anywhere



The V-20 is **battery powered**, so can be **used absolutely anywhere regardless of where the nearest plug socket is.**

## Key features

Lightweight design  
Unit switchable between kg and lb  
Carry handle  
Previous weight reading memory  
Includes Tare and Hold functions


## Built to last

The V-20 is **robustly built** and **designed to last.** This scale like all our veterinary scales, is backed by a **full one year warranty.**

## Specification

Capacity: 20kg  
Graduations: 10g  
Power: Mains/1 x 9V battery  
Base dimensions: 620mm x 350mm x 20mm  
Weight of scale: 2.5kg

# Marsden V-20: Key features at-a-glance



# MARSDEN

Tel: 01709 364296 / 0800 169 2775

Fax: 01709 364293

E-mail: [sales@marsdengroup.co.uk](mailto:sales@marsdengroup.co.uk)

[www.marsden-weighing.co.uk](http://www.marsden-weighing.co.uk)

 MarsdenWeighing

 @MarsdenWeighing

 [www.youtube.com/MarsdenWeighing](http://www.youtube.com/MarsdenWeighing)

 <http://pinterest.com/marsdenweighing>