

MARSDEN

V-22

Marsden Veterinary Scale



Lightweight, portable pet scale for home or surgery use

Highly accurate to 5g

Wipe clean surface

Perfect for small pets

www.marsden-weighing.co.uk



Marsden V-22: A versatile animal scale for home and surgery use

The **Marsden V-22** is a **lightweight animal scale** suitable for a variety of pets. With **simple but essential features**, it's **very easy to use**.

Keep regular check of your pets weight

The V-22 is a great **low cost weighing scale** which is ideal for use in your own home or in veterinary surgeries. This scale **can be used for weighing a vast range of animals from hamsters, cats and dogs**.

The **large LCD display** clearly **shows the weight reading to the nearest 5g** ensuring that the animals are always measured to the highest degree of accuracy.

Due to the compact shape and size of this scale it can be stored absolutely anywhere – like a cupboard or under a bed!

Quick and easy to use

The V-22 is **powered by 4 x AA batteries**. To start using the scale simply press the power button and you are good to go.

It has a **Tare function** for removing unwanted weight like a blanket and a **Memory Recall** button, which allows the user to recall the previous weight reading. The **units button switches the scale from weighing in kg, to weighing in lb**.



Removable tray reveals a platform for smaller pets

The wipe-clean tray is perfect for weighing cats and dogs, but remove it and you have a **smaller weighing platform that's perfect for smaller pets** – like hamsters!

Key features

Wipe-clean surface

Weighing tray reveals smaller platform beneath

Weighs in kg or lb

Previous weight reading memory

Includes Tare and Hold functions

Built to last

The Marsden V-22 is robustly built and designed to last. Like all our veterinary scales, this scale is backed by a **full one year warranty**.

Specification

Capacity: 20kg

Graduations: 5g

Power: 4x AA batteries

Dimensions: 585mm x 335mm x 115mm

Weight of scale: 2.7kg

Marsden V-22: Key features at-a-glance



MARSDEN

Tel: 01709 364296 / 0800 169 2775

Fax: 01709 364293

E-mail: sales@marsdengroup.co.uk

www.marsden-weighing.co.uk

[f](#) MarsdenWeighing

[t](#) @MarsdenWeighing

[You Tube](#) www.youtube.com/MarsdenWeighing

[p](#) <http://pinterest.com/marsdenweighing>