

MARSDEN

V-23

Marsden Portable Pet Scale



Lightweight, portable veterinary scale

25kg capacity; accurate to 10g

Non-slip, wipe clean surface

Perfect for small pets

www.marsden-weighing.co.uk



Marsden V-23: A portable weighing scale for small pets at home or in the surgery

The **V-23** is a **small, lightweight animal scale**, perfect for keeping regular track of a pet's weight. It's been designed as a weighing scale that can be used anywhere, anytime, and safely stored away when not in use.

Keep regular check of your pets weight

The V-23 is a veterinary scale designed for the smaller pet. It can weigh pets up to a maximum of 25kg, with accuracy to the nearest 10g.

Plus, the V-23 has **Tare and auto power off functions**, a non-slip, wipe clean surface, and an easy-to-use **Hold feature**, meaning the weight reading Holds however much your pet moves whilst on the scale.

Take it anywhere you need

The V-23 is battery powered scale, meaning it can be used absolutely anywhere - with **four AAA batteries** providing power for at least 40 hours.

It has a **built-in carry handle** so it can be carried around - and even hung from the wall when not in use.

This wipe clean scale is of very light construction - making it extremely portable. Plus, its small size means it can be stored away when not in use.



Keep a track of changes in your pet's weight

The V-23 stores the previous weight reading and shows the weight difference. This makes it quick and easy to see if your pet's weight has gone up, or down, from their previous weight check.



Key features

- Lightweight design
- Carry handle
- Weighs in kg or lb
- Previous weight reading memory
- Includes Tare and Hold functions



Built to last

The Marsden V-23 is robustly built and designed to last. Like all our veterinary scales, this scale is backed by a **full one year warranty**.

Specification

- Capacity: 25kg
- Graduations: 10g
- Power: 4x AAA batteries
- Dimensions: 320mm x 550mm x 83mm
- Weight of scale: 1.5kg

Marsden V-23: Key features at-a-glance



MARSDEN

Tel: 01709 364296 / 0800 169 2775

Fax: 01709 364293

E-mail: sales@marsdengroup.co.uk

www.marsden-weighing.co.uk

 MarsdenWeighing

 @MarsdenWeighing

 www.youtube.com/MarsdenWeighing

 <http://pinterest.com/marsdenweighing>