

# MARSDEN

## V-250

Large Veterinary Scale



High capacity, super-size veterinary scale

Stainless steel casing for easy cleaning

Low profile 50mm platform

Optional non-slip rubber mat for animal comfort

[www.marsden-weighing.co.uk](http://www.marsden-weighing.co.uk)



## Marsden V-250: A large, stainless steel wash-down weighing platform

The **V-250** is a superbly designed scale with a super-size stainless steel base. This scale is suitable to weigh anything from domestic cats to the largest of dogs!

### Suitable for large goods

A larger veterinary weighing scale without the premium price - that's the V-250. With its large, low profile base, the **V-250 is capable of weighing animals upto a capacity of 250kg.**

The large base is easy for pets to step onto and a rubber mat is also available as an accessory.

### Easy to use functions

The V-250 has been **designed for simple, accurate animal weighing.** There is a **large, clear LCD display** for the weight reading, and four function buttons, For **Tare** (to **remove unwanted weight** from the display) and **Hold** (to **stabilise the weight reading** if the animal is moving on the scale).



## Perfect for surgery waiting areas



The size and simplicity of the V-250 means it's **perfect for veterinary surgery waiting areas**. Plus, **battery-powered operation** means you can position the scale wherever you need it. The V-250 can **also be powered from the mains**.

## Key features

- Stainless steel wipedown base
- Large easy to read LCD display
- Tare and Hold functions
- Weighs in kg or lb
- Wall bracket for mounting indicator



## Get peace of mind with a service contract

Like all Marsden veterinary scales, you can keep your V-250 **accurate and reliable** for longer with a **service contract**. For more information about service contracts please contact us.

## Specification

- Capacity: 250kg
- Graduations: 100g
- Power: Mains or 6 x AA Batteries
- Dimensions: 1090mm x 505mm x 70mm
- Weight of scale: 18kg

# Marsden V-250: Key features at-a-glance



# MARSDEN

Tel: 01709 364296 / 0800 169 2775

Fax: 01709 364293

E-mail: [sales@marsdengroup.co.uk](mailto:sales@marsdengroup.co.uk)

[www.marsden-weighing.co.uk](http://www.marsden-weighing.co.uk)

 MarsdenWeighing

 @MarsdenWeighing

 [www.youtube.com/MarsdenWeighing](http://www.youtube.com/MarsdenWeighing)

 <http://pinterest.com/marsdenweighing>