

MARSDEN

V-24

Robust Marsden Veterinary Scale



Durable pet scale for home or surgery use

Highly accurate to 10g

Hold, Tare and Unit functions

Ultra-lightweight; perfect for small pets

www.marsden-weighing.co.uk



Marsden V-24: An ultra-lightweight, durable veterinary scale for home and surgery use

The **Marsden V-24** is a highly robust, portable animal scale, perfect for weighing small pets. At just 1.55kg, it is **one of the lightest vet scales in our range**. With **simple but essential features**, it's **very easy to use**.

Keep regular check of your pets weight

The V-24 is a **highly accurate veterinary weighing scale** which is ideal for use in your own home or in veterinary surgeries. This scale **can be used for weighing a range of animals from hamsters to cats and dogs**.

The **bright LCD display** clearly **shows the weight reading to the nearest 10g** ensuring that animals are always weighed with a very high degree of accuracy.

The lightweight design means the scale is easy to move and can be stored absolutely anywhere. **Built in overload protection** helps protect the scale from damage.

Quick and easy to use

The V-24 is **powered by 3 x AAA batteries**. To start using the scale simply press the power button and you are good to go. With power for **60 hours' continuous use**, it has the longest battery life of any Marsden veterinary scale.

It has a **Tare function** for removing unwanted weight like a blanket and the **Hold function** maintains the reading on the LCD display even after the pet has left the scale. The **units button switches the scale from weighing in kg, to weighing in lb**.



Removable tray reveals a platform for smaller pets

The V-24 features a **modern, sleek design**. The wipe-clean sculpted tray is perfect for weighing cats and dogs, but remove it and you have a **smaller weighing platform that's perfect for smaller pets** – like hamsters!

Key features

- Wipe-clean surface
- Weighing tray reveals smaller platform beneath
- Weighs in kg or lb
- Bright LCD display
- Includes Tare and Hold functions
- Overload protection
- Portable
- Lightweight

Built to last

The Marsden V-24 is a highly robust vet scale and is designed to last. The scale has an **overload capacity of 100kg**, meaning there is little or no chance of damage from overloading the scale.

Like all our veterinary scales, this scale is backed by a **full one year warranty**.

Specification

- Capacity: 20kg
- Graduations: 10g
- Power: 3x AAA batteries
- Dimensions: 609mm x 326mm x 91mm
- Weight of scale: 1.55kg

Marsden V-24: Key features at-a-glance



MARSDEN

Tel: 01709 364296 / 0800 169 2775

Fax: 01709 364293

E-mail: sales@marsdengroup.co.uk

www.marsden-weighing.co.uk

 MarsdenWeighing

 @MarsdenWeighing

 www.youtube.com/MarsdenWeighing

 <http://pinterest.com/marsdenweighing>