


Weigh time critical and immobile patients instantly, accurately

- The Patient Transfer Scale is a transfer board with an in-built weighing scale. It provides a weight reading whilst transferring a patient from trolley to bed or bed to bed
- A more dignified way of weighing, reducing stress for patients
- It's part of an existing, embedded process: quick and easy, saving time
- Provides an accurate weight when it's critical for safe medication management and monitoring of patients



Introducing the Patient Transfer Scale

MARSDEN



UK PATENT

Patent number 2543835



US PATENT PENDING



EU PATENT PENDING



CANADA PATENT PENDING



HONG KONG PATENT PENDING

The Patient Transfer Scale:

A revolutionary solution for weighing time critical or immobile patients quickly



Key Features

- 93/42/EEC Medical Devices Directive Approved
- 2014/31/EU Non-Automatic Weighing Instrument Directive
- Class III Approved
- Large weight display
- Portable
- Wipe clean surface
- Hold-of-weight feature

Specifications

- Capacity: 250kg
- Graduations: 500g
- Power: Internal rechargeable battery
- Battery life: 50 hours
- Weight: 11kg
- Dimensions: 1800mm x 700mm x 30mm
- Product code: M-999

- The idea of the Patient Transfer Scale (PTS) was born of frustration: An NHS Lanarkshire nurse working in Emergency Care was determined to find an easier way to weigh immobile patients
- Healthcare professionals have a duty of care to ensure a safe environment for patients, whilst maintaining their dignity. The PTS provides clinicians with confidence, providing an instant, accurate weight
- The Patient Transfer Scale means you no longer need to estimate the weight of a patient
- A patient's weight can be critical to delivering the right care and treatment at the right time
- Tested in over 30 hospitals, the PTS' final design is based on feedback from healthcare professionals
- The PTS is lightweight and portable and can be taken to exactly where your patient is

Gillian Taylor

Clinical Consultant/Inventor of PTS
 +44 (0) 7397 787780
 gillian@marsdengroup.co.uk
 @gilliantaym999

Paul McCaig

Sales Consultant
 +44 (0) 7885 998511
 paulm@marsdengroup.co.uk

Marsden Head Office Sales

+44 (0)1709 364296
 www.patienttransferscale.co.uk
 www.marsden-weighing.co.uk

