

MARSDEN

M-125

ITEM NO	DESCRIPTION	PART NO	PRICE £
1	I-300 MAIN BOARD WITH BUTTONS	SE90008001	125.00
2	I-300 LCD DISPLAY	ME71000572	18.00
3	I-300 BATTERY HOLDER	MP00080-23	9.00
4	I-300 BATTERY DOOR	MP04000541	6.00
5	I-300 OVERLAY	MP10004951	22.00
6	I-300 LENS	MP10004951	9.00
7	DP-3810 ADAPTOR 12V 1A (TIP+)	UE091105DGLV01-P	36.00
8	AA BATTERY X 6 (NON RECHARGEABLE)		3.00 EACH
9	TOP HOUSING	MP01001421	30.00
10	BOTTOM HOUSING	MP02000981	31.00
11	BATTERY CONNECTOR	SE00301121	18.00
12	BUTTONS	ME73000831	3.00

ITEM NO	DESCRIPTION	PART NO	PRICE £
1	BLACK PLASTIC HEIGHT MEASURE BRACKET	M125-2501	25.00
2	LOAD CELL ZEMIC L6E3 300KG	L6E3 300KG	75.00
3	FEET 8MM	M125-2502	6.00
4	LOAD CELL LEAD	M125-2503	30.00
5	HEIGHT MEASURE CURSOR	M125-2504	18.00
6	SIDE MOUNT BRACKET	M1252505	15.00

Sales: 0845 130 7330
sales@marsdengroup.co.uk

Service: 01709 364 296
service@marsdengroup.co.uk

Website:
marsden-weighing.co.uk