

MARSDEN

M-430

ITEM NO	DESCRIPTION	PART NO	PRICE £
1	OVERLAY	MP10003751	14.00
2	LCD DISPLAY BOARD	PF06700021	35.00
3	LOAD CELL L6D8(150KG)	ME93000511	60.00
4	A/D BOARD	PF06700011	85.00
5	FEET	MP00600192	5.00
6	BATTERY BOX	MP00001691	8.00
7	BATTERY COVER	AE-0162	6.00
8	BABY TRAY	ME9311115	40.00
9	SLIDING SWITCH ON/OFF	AE-0992	4.00
10	AC JACK SOCKET	ME04000291	7.00
11	ADAPTOR 12V	LGSPSB120100BS	38.00
12	TILE	M400 TILE WHITE	22.00
13	LENS COVER	MP10100181	9.00
15	VR BOARD ASSEMBLY	90016001475	12.00

Sales: 0845 130 7330
sales@marsdengroup.co.uk

Service: 01709 364 296
service@marsdengroup.co.uk

Website:
marsden-weighing.co.uk