

MARSDEN


BP2 Blood Pressure Monitor with ECG



- Blood pressure monitor with integrated ECG
- Bluetooth connectivity for iOS/Android
- Clinically validated device
- Easy to use with only two buttons
- Real time measurements
- Portable and on-the-go

Marsden BP2:

Blood Pressure Monitor with ECG



- Blood pressure monitor with integrated ECG: help you monitor your heart health in your own home
- Bluetooth connectivity: By connecting with our free app, you can monitor your results over time to look for trends and irregularities, such as hypertension and arrhythmia, and help you to discuss your concerns with a medical professional. Results can also be shared from the app to a medical professional or a family member
- Clinically validated: CE Marked and FDA approved, so you can be confident that the results you take and maybe share with medical professionals are consistent and reliable.
- Real time measurements: Take your blood pressure and ECG and see them in real time in the application if desired.
- Easy to use: only 2 buttons to use, follow the on-screen instructions, check your health status within 30 seconds. Instant results are displayed on the device's OLED screen. The built-in memory can store 10 sets of ECG waveforms and 50 sets of blood pressure data.
- Portable and on-the-go: the compact, one-piece design makes it easy to monitor your blood pressure in the comfort of your home, while traveling or at work.

Weight: 200 g

Warranty Period: 1 Year

Power Supply: USB Charger

Battery Life: 50 hours of continuous use