

MARSDEN

H-628 Free Standing Height Measure



- Measuring range: 75cm - 215cm
- Graduations of 1mm
- Aluminium height rod for hospitals
- Free standing with base

Marsden H-628: Free Standing Height Measure

- Measuring range: 75cm - 215cm
- Graduations of 1mm
- Aluminium height rod for hospitals
- Free standing with base
- Easy to set up and use
- Suitable for GP surgeries
- Sliding cursor
- 1 year warranty as standard



Free standing design

The H-628 is a free standing height measure, which does not need to be used in conjunction to a weighing scale, or mounted on a wall in order to gain an accurate height.

We designed this height measure after receiving feedback from our network of hospitals and GP surgeries asking for something which was easy to operate and could be transported from room to room.

Ideal for GP surgeries

During the design stages of this height measure, we found that GP surgeries had an increasing requirement for a height measure to calculate BMI.

The device also had to be a space saving design which would not take up too much room in small surgeries and with a lightweight construction to easily move around and store.

Easy to use

The H-628 is a very simple and easy to use height measure. From the box, the height measure is easy to construct, taking only a few minutes.

Once assembled, the patient needs to simply stand on the base and the healthcare professional will carefully lower the cursor until it is sitting level on top of the patients head.

Measuring Range: 75cm - 215cm

Graduations: 1mm

Base Dimensions: 350mm x 380mm x 50mm

Weight: 0.6kg