

MARSDEN

HM-200D Digital Height Measure



- Measuring range: 120cm - 200cm
- Graduations: 1mm
- Telescopic height measure made of aluminium
- Fits to Marsden column scales

Marsden HM-200D: Digital Height Measure

- Measuring range: 120cm - 200cm
- Graduations: 1mm
- Telescopic height measure made of aluminium
- Fits to Marsden column scales
- Digital display to show the height
- Incorporates folding cursor for safety
- Patients height read at eye level
- Covered by 1 year warranty



Digital display for accuracy

The HM-200D is a height measure with a small digital display which shows you a clear and accurate height once you have pulled the cursor down.

This height measure has proven to be incredibly popular as it removes any guesswork as to where the indicator is situated and it is easy to see from any height thanks to the LED display.

Attach to our column scale

As well as the ability to attach this height measure to a wall, it can also be mounted to our popular M-110 and M-125 column weighing scales.

If the height measure is purchased at the same time as your column scale, the device will be delivered as a complete solution. Alternatively, you can add the HM-200D as an addition to an existing device.

Easy to use

The HM-200D is an incredibly simple and easy to use height measure. This device was designed with making the measurement process as easy for the user as possible.

Whilst your patient is stood on the scale, simply turn on the power button on the height measure, then lower the cursor until it is sat on the head of your patient. The reading will now clearly be displayed.

Measuring Range: 120cm - 200cm

Graduations: 1mm

Dimensions: 1000mm x 75mm x 90mm

Weight: 0.6kg