

MARSDEN

HM-201M Manual Height Measure



- Measuring range: 60cm - 200cm
- Graduations: 1mm
- Telescopic height measure made of aluminium
- Fits to Marsden column scales

Marsden HM-201M: Manual Height Measure

- Measuring range: 60cm - 200cm
- Graduations: 1mm
- Telescopic height measure made of aluminium
- Fits to Marsden column scales
- Can be fitted to a wall for saving space
- Incorporates folding cursor for safety
- Patients height read at eye level
- Covered by 1 year warranty



Mount on a wall

The HM-201M is a multi-functional height measure which can be easily attached to a wall in order to save space in hospitals and GP surgeries.

The height measure will come as standard with all the fixtures and fittings which you may need in order to attach to the wall and the device can then be easily used and is a space saving device.

Attach to our column scale

As well as the ability to attach this height measure to a wall, it can also be mounted to our popular M-110 and M-125 column weighing scales.

If the height measure is purchased at the same time as your column scale, the device will be delivered as a complete solution. Alternatively, you can add the HM-201M as an addition to an existing device.

Easy to use

Whether the HM-201M height measure is attached to the wall or a column scale, the device is incredibly simple to use and easy to operate.

Simply ask your patient to stand near the height measure and then once they are in position, the healthcare worker can carefully and slowly lower the cursor into position until it is sitting on top of their head.

Measuring Range: 60cm - 200cm

Graduations: 1mm

Dimensions: 1000mm x 90mm x 55mm

Weight: 0.6kg