


MARSDEN

HM-210D Digital Wall Mounted Height Measure



- Measuring range: 100cm - 210cm
- Graduations: 1mm
- Wall mounted height measure
- Accurate digital display

Marsden HM-210D: Wall Mounted Digital Height Measure

- Measuring range: 100cm - 210cm
- Graduations: 1mm
- Fits to a wall easily
- Digital display for accurate readings
- Simple and easy to use
- Incorporates folding cursor for safety
- Patients height read at eye level
- Covered by 1 year warranty



Digital display for accuracy

The HM-210D is a height measure with a small digital display which shows you a clear and accurate height once you have pulled the cursor down.

This height measure has proven to be incredibly popular as it removes any guesswork as to where the indicator is situated and it is easy to see from any height thanks to the LED display.

Multi-use height measure

The HM-210D is a wall mounted height measure which was designed to weigh a broad range of individuals including children, adults and the elderly.

Simply mount the height measure on the wall of your clinic or GP surgery using the two hinges which are supplied as standard, as you are ready to begin getting accurate readings for your patients.

Easy to use

The HM-210D is a simple and easy to use height measure, which takes seconds to use, once it is in position and affixed to the wall where it will be situated.

A patient needs to simply walk up to the height measure, removing shoes if necessary and the healthcare worker will then slowly move the cursor down until it sits on the patients head.

Measuring Range: 100cm - 210cm

Graduations: 1mm

Dimensions:

2140mm x 70mm x 100mm

Weight:

0.6kg