MARSDEN

HM-80P Baby/Toddler Measuring Rod



- · Lightweight plastic construction
- Sliding cursor for accurate measuring
- 1mm increments
- Suitable for babies and toddlers

Marsden HM-80P:

Height Measure for Babies

- Class III Approved
- High level of accuracy
- Perfect for use in neo-natal units
- Mains operated or AA batteries with over 600 hours of use
- Auto power-off function helps to save battery power
- One-piece baby scale to cradle the baby

- Lightweight design making it ideal for health visitors and midwives
- Large, bright display which is easy to read
- Adjustable levelling feet to ensure that the scale is stable
- Hold and Tare functions to enable accurate and easy weighing
- Optional rucksack and carry case available for transportation
- Covered by the Marsden 4 year warranty



Multi-use Height Measure

The Marsden HM-80P is suitable for both babies and toddlers thanks to the simplistic design and easy to use sliders to adjust to the length of the infant.

With a max measuring range of 80cm, this height measure can accommodate a new-born up to 18-24 months of age, with a precise accuracy of only 1mm. Ideal for close monitoring.

Extreme accuracy

With 1mm increments, height stopper and red arrows for precision, the Marsden HM-80P is one of the most accurate height measures available on the market.

Whatever the heigh measurement is needed for, you can guarantee that it will be precise and accurate with the HM-80P. Highly popular with medical professionals and at home.

Easy to use

To use the height measure, simply lay the baby or toddler down on a safe, flat surface and position the HM-80P height measure around the infant.

Once the height measure and infant are both in the correct position, place one of the cursors level with the baby's head, then adjust the other cursor until it is level with their feet, their height will then be shown.

Measuring Range:	10cm - 80cm	Dimensions:	880mm x 100mm x 330mm
Graduations:	1mm	Weight:	0.6kg