

MARSDEN

M-100 Column Scale with Height Measure



- Class III Approved
- Accurate to 50g
- Tough, robust and hard wearing base
- Corner design for saving space
- Suitable for GPs and hospitals
- Optional Wi-Fi or Bluetooth available

Marsden M-100:

Column Scale with Integrated Height Measure

- Class III Approved
- Accurate to 50g
- Perfect for use in the GP surgery or hospital
- Powered by rechargeable battery pack
- 55 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Large, robust base for weighing bariatric patients
- Easy to read display
- Hold and Tare functions to enable accurate and easy weighing
- Optional Bluetooth and Wi-Fi available
- Covered by the Marsden 4 year warranty



Simple weighing

The M-100 is a simple and easy to use column scale, constructed to the highest quality to ensure that it is suitable for weighing patients up to 300kg.

The scale has been designed to be hard wearing making it ideal for heavy use in busy traffic areas, such as general wards, outpatient departments and busy GP surgeries.

Optional connectivity options

The M-100 is available with either Bluetooth or Wi-Fi connectivity to ensure that our customers are able to transfer data seamlessly or input into EPR systems.

We would simply need to insert a connectivity module at the time of order and then provide you with our standard protocols for easy and simple integration into your existing system.

Easy to use

The M-100 features from our best-selling DP-3810 indicator. This device features from a full numeric keypad for easy BMI calculation and simple function keys.

The On/Off button ensures that you have full control over the device and you can prolong the battery life. Hold enables you to gain readings for unsteady patients and Tare allows you to deduct unwanted weight.

Capacity:	300kg	Power Supply:	Mains/Rechargeable
Graduations:	50g<150kg>100g	Battery Life:	55 hours of continuous use
Functions:	Hold, Portable, Tare	Approvals:	Class III Approved
Dimensions:	460mm x 440mm x 1270mm	Optional Extras:	Bluetooth/Wi-Fi, Digital Height Measure, Printer
Weight of Scale:	15kg	Warranty Period:	4 Years