

MARSDEN

M-125 Entry Level Column Scale



- Class III and MDD Approved
- Accurate to 100g
- Tough, robust and hard wearing base
- Simple and easy to use
- Suitable for GPs and hospitals
- BMI and BSA function

Marsden M-125:

Entry Level Column Scale

- Class III and MDD Approved
- Accurate to 100g
- Perfect for use in the GP surgery or hospital
- Powered by 6 x AA batteries
- 180 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Optional carry case available to aid portability
- Easy to read remote display
- Hold and Tare functions to enable accurate and easy weighing
- Made from reinforced aluminium structure
- Covered by the Marsden 8 year warranty



Simple weighing

The M-125 is a simple and easy to use stand on scale, made from reinforced aluminium, which is ideal for weighing all patients up to a capacity of 250kg.

The scale has been designed to be hard wearing making it ideal for heavy use in busy traffic areas, such as general wards, outpatient departments and busy GP surgeries.

With BSA function

BSA function calculates the surface area of the human body. BSA is commonly used by NHS DEXA scanners to help calculate body composition and medication dosages.

We believe Marsden is the first digital weighing scales manufacturer in the world to include BSA calculation function on a weighing scale. More can be found about BSA on our website.

Easy to use

The M-125 features our popular I-300 indicator, which is very simple to use due to the 5 buttons and a bright, clear to read display.

The On/Off button ensures that you have full control over the device. Hold key enables you to gain weights even for the patient is not steady on the scale and the tare key allows you to deduct any unwanted weight.

Capacity:	250kg	Power Supply:	6 x AA Batteries
Graduations:	100g	Battery Life:	180 hours of continuous use
Functions:	Hold, Portable, Tare	Approvals:	Class III and MDD Approved
Dimensions:	440mm x 310mm x 940mm	Optional Extras:	Manual Height Measure
Weight of Scale:	9kg	Warranty Period:	8 Years