

# MARSDEN

## M-150 Gym and Fitness Scale



- Capacity of 300kg
- Accurate to 200g
- Height measure included as standard
- Weighs in stones or kilos
- BMI calculation
- Suitable for gyms and health clubs

# Marsden M-150:

## Gym and Fitness Column Scale

- Capacity of 300kg/47st
- Accurate to 200g
- Perfect for use in gyms and health clubs
- Powered by rechargeable battery
- Simple and easy to use
- Optional carry case available to aid portability
- Hold and Tare functions to enable accurate and easy weighing
- Bright, clear LED display
- Covered by the Marsden 1 year warranty
- Unit switchable



### Simple weighing

The M-150 is a simple and easy to use stand on scale. This scale is non-approved with switchable units, making it perfect for use in gyms and health clubs.

The scale has a capacity of 300kg or 47 stone and is accurate to graduations of 200g, which makes it suitable for customers of all shapes and sizes to track their weightloss journey.

### With BMI calculation

The M-150 is an all-in-one solution for people wanting to monitor their weight as it comes as standard with BMI function, which is a universal indicator for health.

Use the height measure which comes as standard and then input that figure into the M-150 using the arrow keys. The scale will then automatically generate the BMI and it will appear on screen.

### Easy to use

The M-150 has proven to be a highly popular gym and fitness scale, thanks to the four simple to use keys and large easy to read digital display.

The On/Off button ensures that you have full control over the device. Hold enables you to gain weights even for the patient is not steady on the scale and the tare key allows you to deduct any unwanted weight.

<b>Capacity:</b>	300kg/47st	<b>Power Supply:</b>	Rechargeable Battery
<b>Graduations:</b>	200g	<b>Battery Life:</b>	40 hours of continuous use
<b>Functions:</b>	Hold, Portable, Tare, BMI	<b>Approvals:</b>	Non-Approved
<b>Dimensions:</b>	970mm x 310mm x 540mm	<b>Optional Extras:</b>	N/A
<b>Weight of Scale:</b>	9.2kg	<b>Warranty Period:</b>	1 Year