

MARSDEN

M-200 High Capacity Chair Scale



- Class III and MDD Approved
- Accurate to 50g
- Tough, robust and hard wearing scale
- Full numeric keypad for BMI calculation
- Suitable for Hospitals
- Optional Wi-Fi or Bluetooth available

Marsden M-200:

Professional High Capacity Chair Scale

- Class III and MDD Approved
- Accurate to 50g
- Perfect for use in hospitals
- Powered by rechargeable battery pack
- 55 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Hard-wearing scale, built to last
- Easy to read display
- Folding hand rests and footrests
- Optional Bluetooth and Wi-Fi available
- Covered by the Marsden 4 year warranty



Simple weighing

The M-200 is a simple and easy to use chair scale, which has been manufactured from the highest quality materials to ensure longevity and accuracy.

This scale is suitable for use in many hospital departments, including outpatients departments, general wards and is suitable for heavy use and high traffic areas with many patients.

Extra wide seat

In order to help accommodate more patients, we designed the M-200 to have a larger seat than standard chair scales to provide more room for patients.

The M-200 is also available in capacities of 250kg and 300kg, to enable us to provide our customers with highly accurate readings for their patients, making monitoring, diagnosis and treatment easier.

Optional connectivity options

The M-200 is available with either Bluetooth or Wi-Fi connectivity to ensure that our customers are able to transfer data seamlessly or input into EPR systems.

We would simply need to insert a connectivity module at the time of order and then provide you with our standard protocols for easy and simple integration into your existing system.

Capacity:	250kg or 300kg	Power Supply:	Mains/Rechargeable
Graduations:	50g<150kg>100g	Battery Life:	55 hours of continuous use
Functions:	Hold, Portable, Tare, BMI	Approvals:	Class III and MDD Approved
Dimensions:	1000mm x 690mm x 990mm	Optional Extras:	Printer
Weight of Scale:	27kg	Warranty Period:	4 Years