

MARSDEN

M-210 Professional Chair Scale



- Class III Approved
- Accurate to 50g
- Tough, robust and hard wearing scale
- Full numeric keypad for BMI calculation
- Suitable for hospitals and care homes
- Optional Wi-Fi or Bluetooth available

Marsden M-210:

Professional Chair Scale

- Class III Approved
- Accurate to 50g
- Perfect for use in hospitals and care homes
- Powered by rechargeable battery pack
- 40 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Hard-wearing scale, built to last
- Easy to read display
- Adjustable armrests and footrests
- Optional Bluetooth and Wi-Fi available
- Covered by the Marsden 8 year warranty



Simple weighing

The M-210 is a simple and easy to use chair scale, which has been manufactured from the highest quality materials to ensure longevity and accuracy.

This scale is suitable for use in many hospital departments, including outpatients departments, general wards and is suitable for heavy use and high traffic areas with many patients.

Folding armrests/footrests

Manoeuvring patients is different for every individual. For this reason, the M-210 to have adjustable armrests and footrest, to make the process as easy as possible.

This chair scale has gone through rigorous testing by various hospitals and care homes throughout the country and with their feedback, we have ensured that the scale is not only functional, but also comfortable and safe for patients.

Optional connectivity options

The M-210 is available with either Bluetooth or Wi-Fi connectivity to ensure that our customers are able to transfer data seamlessly or input into EPR systems.

We would simply need to insert a connectivity module at the time of order and then provide you with our standard protocols for easy and simple integration into your existing system.

Capacity:	250kg	Power Supply:	Mains/Rechargeable
Graduations:	50g<150kg>100g	Battery Life:	40 hours of continuous use
Functions:	Hold, Portable,, BMI	Approvals:	Class III Approved
Dimensions:	880mm x 600mm x 1010mm	Optional Extras:	Printer
Weight of Scale:	18kg	Warranty Period:	8 Years