

# MARSDEN

## M-225 Easy To Use Chair Scale



- Class III Approved
- Accurate to 100g
- Large, clear easy to read display
- Simple and easy to use
- Suitable for hospitals and care homes
- BMI and BSA function

# Marsden M-225:

## Entry Level Chair Scale

- Class III Approved
- Accurate to 100g
- Perfect for use in hospitals and care homes
- Powered by 6 x AA batteries
- 180 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Easy to read remote display
- Hold and Tare functions to enable accurate and easy weighing
- USB Type B port in indicator
- Large, 30mm weight display
- Covered by the Marsden 4 year warranty



### Simple weighing

The M-225 is a simple and easy to use chair scale, made from reinforced aluminium, which is ideal for weighing all patients up to a capacity of 250kg.

The scale has been designed to be hard wearing making it ideal for heavy use in busy traffic areas, such as general wards, outpatient departments and busy GP surgeries.

### With BSA function

BSA function calculates the surface area of the human body. BSA is commonly used by NHS DEXA scanners to help calculate body composition and medication dosages.

We believe Marsden is the first digital weighing scales manufacturer in the world to include BSA calculation function on a weighing scale. More can be found about BSA on our website.

### Easy to use

The M-225 features our popular I-300 indicator, which is very simple to use due to the 5 simple buttons and the large, clear to read display..

The On/Off button ensures that you have full control over the device. Hold enables you to gain weights for the patient is not steady on the scale and tare key allows you to deduct any unwanted weight.

<b>Capacity:</b>	250kg	<b>Power Supply:</b>	Mains or 6 x AA Batteries
<b>Graduations:</b>	100g	<b>Battery Life:</b>	180 hours of continuous use
<b>Functions:</b>	Hold, Portable, Tare, BMI, BSA	<b>Approvals:</b>	Class III Approved
<b>Dimensions:</b>	800mm x 600mm x 1010mm	<b>Optional Extras:</b>	N/A
<b>Weight of Scale:</b>	18kg	<b>Warranty Period:</b>	4 Years