

MARSDEN

M-250 Chair Scale with Stand Assist



- Class III and MDD Approved
- Raising seat, to help patients safely to their feet
- Accurate to 100g
- Tough, robust and hard wearing scale
- Full numeric keypad for BMI calculation
- Optional Wi-Fi or Bluetooth available

Marsden M-250:

Professional Chair Scale with Stand Assist

- Class III and MDD Approved
- Raising seat to help ease patients steadily to their feet
- Accurate to 100g
- Perfect for use in hospitals and care homes
- Powered by rechargeable battery pack
- 55 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Hard-wearing scale, built to last
- Easy to read display
- Optional Bluetooth and Wi-Fi available
- Covered by the Marsden 4 year warranty



Powered seat to help patients

Chair scales are essential when a patient or resident needs to be seated whilst being weighed, especially those with limited mobility.

The M-250 changes this. Using the remote control, the seat steadily tilts, slowly raising patients to an upright position, easing them steadily onto their feet with the help of a medical professional.

Simple weighing

The M-250 is a simple and easy to use chair scale, which has been manufactured from the highest quality materials to ensure longevity and accuracy.

This scale is suitable for use in many hospital departments, including outpatients departments, general wards and is suitable for heavy use and high traffic areas with many patients.

Optional connectivity options

The M-250 is available with either Bluetooth or Wi-Fi connectivity to ensure that our customers are able to transfer data seamlessly or input into EPR systems.

We would simply need to insert a connectivity module at the time of order and then provide you with our standard protocols for easy and simple integration into your existing system.

Capacity:	250kg	Power Supply:	Mains/Rechargeable
Graduations:	100g	Battery Life:	55 hours of continuous use
Functions:	Hold, Portable, Tare, BMI	Approvals:	Class III and MDD Approved
Dimensions:	1025mm x 575mm x 990mm	Optional Extras:	Printer
Weight of Scale:	29kg	Warranty Period:	4 Years