

# MARSDEN

## M-300 Lightweight Baby Scale



- Ultra-lightweight Class III Approved scale
- Highly accurate to 2g graduations
- 600 hours' battery life
- Perfect for maternity wards and midwives
- Moulded design for all baby sizes
- Optional carry case available

# Marsden M-300: Ultra-Lightweight Baby Scale

- Class III Approved
- High level of accuracy
- Perfect for use in neo-natal units
- Mains operated or AA batteries with over 600 hours of use
- Auto power-off function helps to save battery power
- One-piece baby scale to cradle the baby
- Lightweight design making it ideal for health visitors and midwives
- Large, bright display which is easy to read
- Adjustable levelling feet to ensure that the scale is stable
- Hold and Tare functions to enable accurate and easy weighing
- Optional rucksack and carry case available for transportation
- Covered by the Marsden 4 year warranty



## Ultimate portability

The M-300 is the most portable baby scale of its kind in the Marsden range. We designed this scale to provide a lightweight solution for midwives and health visitors.

Therefore, the M-300 weighs just 3.1kg and is available with a carry case or rucksack for use in the community. Additionally, thanks to picostrain technology, just two AA batteries will provide over 600 hours' use.

## Extreme accuracy

The M-300 is accurate to just 2g. This means even the smallest changes in a baby's weight are detected, and it's suitable for weighing for critical treatment.

The Tare feature makes it quick and easy to remove the weight of additional items such as blankets or swaddling, so you can get a true reading of the baby's weight, ensuring accurate monitoring and treatment.

## Built to last

The M-300 is a weighing device which is built to last and we back this by providing this scale with our 4 year warranty as standard.

If you would like additional peace of mind, we would recommend a Marsden service contract to ensure that your device is weighing within the recommended guidelines.

<b>Capacity:</b>	15kg	<b>Power Supply:</b>	Internal Rechargeable Battery or 2 x AA Batteries
<b>Graduations:</b>	2g<6kg>5g	<b>Battery Life:</b>	600 hours of continuous use
<b>Functions:</b>	Hold, Portable, Tare	<b>Approvals:</b>	Class III Approved
<b>Dimensions:</b>	370mm x 590mm	<b>Optional Extras:</b>	AC Adaptor, Carry Case, Rucksack
<b>Weight of Scale:</b>	3kg	<b>Warranty Period:</b>	4 Years