

# MARSDEN

## M-310 Handheld Sling Baby Scale



- Class III Approved
- Accurate to 10g
- Ultra lightweight and portable
- Perfect for health visitors and midwives
- Suitable for weighing babies at home or surgery
- Optional carry case available

# Marsden M-310: Handheld Sling Baby Scale

- Class III Approved
- Accurate to 10g
- Perfect for use the community or the surgery
- Powered by 6 x AAA batteries
- 80 hours of continuous use
- Auto power-off function helps to save battery power
- Ultra lightweight and highly portable
- Optional carry case available to aid portability
- Large, bright display which is easy to read
- Sold with reusable sling
- Hold and Tare functions to enable accurate and easy weighing
- Covered by the Marsden 8 year warranty



## Ultra lightweight and portable

Something that we have heard many times over the years since we began supplying the NHS in the early 1970's is that the baby scales are too heavy for community.

For this exact reason, we brought the highly popular M-310 to the market. Weighing only 0.45kg and with an optional carry case available, this is the scale of choice for many midwives and health visitors.

## Ideal for the community

To make this scale suitable for midwives and health visitors, the scale is lightweight and incredibly portable, thanks to the optional carry case (CC-310).

We designed the carry case to ensure that you are able to contain all of the attachments, such as the sling and hook, in one easy to carry holder. Plus, there is a convenient strap on the top.

## Easy to use

With three simple buttons and easy to read display, the M-310 could not be simpler to use and this is why it is one of our best-selling products.

The On/Off button ensures that you have full control over the device. Hold enables you to gain weights even for the most active babies and the tare keys allows you deduct the weight of the sling or any clothing.

<b>Capacity:</b>	10kg	<b>Power Supply:</b>	6 x AAA Batteries
<b>Graduations:</b>	10g	<b>Battery Life:</b>	80 hours of continuous use
<b>Functions:</b>	Hold, Portable, Tare	<b>Approvals:</b>	Class III Approved
<b>Dimensions:</b>	200mm x 80mm x 105mm	<b>Optional Extras:</b>	Carry Case
<b>Weight of Scale:</b>	0.45kg	<b>Warranty Period:</b>	8 Years