

# MARSDEN

## M-400 Portable Baby Scale



- Class III and MDD Approved
- Accurate to 5g
- Lightweight and portable
- Perfect for health visitors and midwives
- Suitable for weighing babies at home or surgery
- Optional carry case available

# Marsden M-400:

## Portable Baby Scale

- Class III and MDD Approved
- Accurate to 5g
- Perfect for use in the community or the surgery
- Powered by 6 x AA batteries
- 100 hours of continuous use
- Auto power-off function helps to save battery power
- Lightweight and highly portable
- Optional carry case available to aid portability
- Easy to read display
- Hold and Tare functions to enable accurate and easy weighing
- Made from reinforced aluminium structure
- Covered by the Marsden 4 year warranty



### Versatile weighing

The M-400 is a versatile baby scale which can be used to weigh both babies and toddlers up to 20kg. Slide off the tray to reveal a platform for toddlers to stand.

The scale is constructed from reinforced aluminium to ensure that the structure is strong and the scale will be reliable and can be used for many different people, for a long duration of time.

### Ideal for the community

We understand that babies and toddlers are weighed at home, just as often as they are in a medical environment and because of this the M-410 is lightweight and portable.

This device does have an optional carry case available (CC-400), to make transporting the scale even easier for community midwives and health visitors. It also adds a layer of protection when in transit.

### Easy to use

With three simple buttons and a large easy to read display, the M-410 could not be simpler to use and this is why it is one of our best-selling products.

The On/Off button ensures that you have full control over the device. Hold enables you to gain weights even for the most active babies and the tare keys allows you deduct the weight of the sling or any clothing.

<b>Capacity:</b>	20kg	<b>Power Supply:</b>	Mains/6 x AA Batteries
<b>Graduations:</b>	5g<10kg>10g	<b>Battery Life:</b>	100 hours of continuous use
<b>Functions:</b>	Hold, Portable, Tare	<b>Approvals:</b>	Class III and MDD Approved
<b>Dimensions:</b>	585mm x 280mm x 120mm	<b>Optional Extras:</b>	Bluetooth/Wi-Fi, Carry Case
<b>Weight of Scale:</b>	4kg	<b>Warranty Period:</b>	4 Years