

MARSDEN

M-410 Portable Baby and Toddler Scale



- Class III and MDD Approved
- Accurate to 10g
- Lightweight and portable
- Perfect for health visitors and midwives
- Suitable for weighing babies at home or surgery
- Optional carry case available

Marsden M-410:

Portable Baby and Toddler Scale

- Class III and MDD Approved
- Accurate to 10g
- Perfect for use in the community or surgery
- Powered by 6 x AA batteries
- 100 hours of continuous use
- Auto power-off function helps to save battery power
- Lightweight and highly portable
- Optional carry case available to aid portability
- Easy to read display
- Hold and Tare functions to enable accurate and easy weighing
- Made from reinforced aluminium structure
- Covered by the Marsden 4 year warranty



Versatile weighing

The M-410 is a versatile baby scale which can be used to weigh babies, toddlers and children up to 50kg. Slide off the tray to reveal a platform for toddlers to stand.

The scale is constructed from reinforced aluminium to ensure that the structure is strong and the scale will be reliable and can be used for many different people, for a long duration of time.

Ideal for the community

We understand that babies and toddlers are weighed at home, just as often as they are in a medical environment and because of this the M-410 is lightweight and portable.

This device does have an optional carry case available (CC-400), to make transporting the scale even easier for community midwives and health visitors. It also adds a layer of protection when in transit.

Easy to use

With three simple buttons and a large easy to read display, the M-410 could not be simpler to use and this is why it is one of our best-selling products.

The On/Off button ensures that you have full control over the device. Hold enables you to gain weights even for the most active babies and the tare keys allows you deduct the weight of the sling or any clothing.

Capacity:	50kg	Power Supply:	Mains/6 x AA Batteries
Graduations:	10g<20kg>20g	Battery Life:	100 hours of continuous use
Functions:	Hold, Portable, Tare	Approvals:	Class III and MDD Approved
Dimensions:	585mm x 280mm x 120mm	Optional Extras:	Bluetooth/Wi-Fi, Carry Case
Weight of Scale:	4kg	Warranty Period:	4 Years