

MARSDEN

M-425 Stand On Scale with Remote Display



- Class III and MDD Approved
- Accurate to 200g
- Lightweight and portable
- Ideal for community health visitors
- Suitable for GPs, hospitals and the home
- Optional carry case available

Marsden M-425:

Digital Stand On Scale with Remote Display

- Class III and MDD Approved
- Accurate to 50g
- Perfect for use in the community, surgery or hospital
- Powered by 6 x AA batteries
- 100 hours of continuous use
- Auto power-off function helps to save battery power
- Lightweight and highly portable
- Optional carry case available to aid portability
- Easy to read display
- Hold and Tare functions to enable accurate and easy weighing
- Made from reinforced aluminium structure
- Covered by the Marsden 4 year warranty



Simple weighing

The M-425 is a heavy duty stand on scale, made from reinforced aluminium, which is ideal for weighing adults as it is accurate to a minimum of 50g.

The scale features from a 'tap-on' function, which means that the user does not need to bend down to operate the scale, and you can simply tap the device with your foot and it will power on.

Remote display as standard

The M-425 comes with a remote display fitted as standard. This means that when you are operating the scale in sensitive environments, the weight reading is discreet.

The remote display also means that the user does not have to bend down and see the weight reading on the display, as it is easy to see on the screen and the user is even able to operate the device.

Easy to use

With three simple buttons, the M-425 could not be simpler to use and this is why it is one of our best-selling products to medical professionals and home users.

The On/Off button ensures that you have full control over the device. Hold enables you to gain weights even for the most active babies and the tare keys allows you deduct the weight of the sling or any clothing.

| | | | |
|-------------------------|----------------------|-------------------------|-----------------------------|
| Capacity: | 220kg | Power Supply: | Mains/6 x AA Batteries |
| Graduations: | 50g<150kg>100g | Battery Life: | 100 hours of continuous use |
| Functions: | Hold, Portable, Tare | Approvals: | Class III and MDD Approved |
| Dimensions: | 320mm x 310mm x 60mm | Optional Extras: | Bluetooth/Wi-Fi, Carry Case |
| Weight of Scale: | 4kg | Warranty Period: | 4 Years |