

# MARSDEN

## M-510 Premium Stand On Scale with BMI



- Class III and MDD Approved
- Accurate to 100g
- Lightweight and portable
- Ideal for outpatient departments
- Suitable for GPs and hospitals
- Optional carry case available

# Marsden M-510:

## Premium Stand On Scale with BMI Function

- Class III and MDD Approved
- Accurate to 100g
- Perfect for use in the GP surgery or hospital
- Powered by rechargeable battery
- 40 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Optional carry case available to aid portability
- Easy to read remote display
- Hold and Tare functions to enable accurate and easy weighing
- Made from reinforced aluminium structure
- Covered by the Marsden 8 year warranty



### Simple weighing

The M-510 is a premium stand on scale, made from reinforced aluminium, which is ideal for weighing adults as it is accurate to a minimum of 100g.

The scale has proven to be a best-seller at Marsden over the years since we first brought this product to market and that is due to the simplicity of the scale and the easy to use functions.

### Remote display as standard

The M-510 comes with a remote display fitted as standard. This means that when you are operating the scale in sensitive environments, the weight reading is discreet for only the user to see.

The remote display also means that the user does not have to bend down and see the weight reading on the display, as it is easy to see on the device.

### Easy to use

The M-510 uses our popular DP-3810 indicator, which features a full numeric keypad, easy to use functions and a clear display.

The On/Off button ensures that you have full control over the device. Hold enables you to gain weights even for the most active babies and the tare keys allows you deduct the weight of the sling or any clothing.

<b>Capacity:</b>	250kg	<b>Power Supply:</b>	Mains/Rechargeable
<b>Graduations:</b>	100g	<b>Battery Life:</b>	40 hours of continuous use
<b>Functions:</b>	Hold, Portable, Tare	<b>Approvals:</b>	Class III and MDD Approved
<b>Dimensions:</b>	360mm x 310mm x 90mm	<b>Optional Extras:</b>	Bluetooth/Wi-Fi, Printer, Carry Case
<b>Weight of Scale:</b>	6kg	<b>Warranty Period:</b>	8 Years