


MARSDEN

M-530 High Capacity Bariatric Scale



- Class III Approved
- Accurate to 50g
- Extra wide base for weighing bariatric patients
- Ideal for outpatient departments
- Suitable for GPs and hospitals
- Optional carry case available

Marsden M-530:

High Capacity Bariatric Scale

- Class III Approved
- Accurate to 50g
- Perfect for use in the GP surgery or hospital
- Powered by rechargeable battery
- 55 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Optional carry case available to aid portability
- Easy to read remote display
- Hold and Tare functions to enable accurate and easy weighing
- Made from reinforced aluminium structure
- Covered by the Marsden 4 year warranty



Simple weighing

The M-530 is a high capacity stand on scale, made from reinforced aluminium, which is ideal for weighing bariatric patients due to the 300kg or 500kg capacity.

The scale has been designed to be ultra-low profile to ensure that it is easy for patients of all ages to carefully step onto, making it suitable for use with geriatric patients in a hospital environment.

Remote display as standard

The M-530 comes with a remote display fitted as standard. This means that when you are operating the scale in sensitive environments, the weight reading is discreet.

The remote display also means that the user does not have to bend down and see the weight reading on the display, as it is easy to see on the device. Plus, the device is entirely controlled via the remote display.

Easy to use

The M-530 uses our popular DP-3810 indicator, which features a full numeric keypad, easy to use functions and a clear display.

The On/Off button ensures that you have full control over the device. Hold enables you to gain weights even for the patient is not steady on the scale and the tare key allows you to deduct any unwanted weight.

Capacity:	300kg or 500kg	Power Supply:	Mains/Rechargeable
Graduations:	50g<150kg>100g / 100g<300kg>200g	Battery Life:	55 hours of continuous use
Functions:	Hold, Portable, Tare	Approvals:	Class III Approved
Dimensions:	595mm x 400mm x 50mm	Optional Extras:	Bluetooth/Wi-Fi, Carry Case, Printer
Weight of Scale:	8kg	Warranty Period:	4 Years