

MARSDEN

M-550BT Highly Accurate Stand On Scale



- Class III and MDD Approved
- Accurate to 200g
- Wide base to accommodate all patients
- Tap-on function
- Suitable for GP surgeries
- Bluetooth option available for connectivity

Marsden M-550:

Entry Level, Highly Accurate Floor Scale

- Class III and MDD Approved
- Accurate to 200g
- Perfect for use in the GP surgery
- Powered by 4 x batteries
- 100 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Optional carry case available to aid portability
- Easy to read remote display
- Hold and Tare functions to enable accurate and easy weighing
- Optional Bluetooth available to enable with data transfer
- Covered by the Marsden 4 year warranty



Simple weighing

The M-550 is a simple and easy to use stand on scale. Simply tap the scale on, using your foot and the scale is then ready to weigh your patient.

The scale has been designed to be hard wearing making it ideal for heavy use in busy traffic areas, but should be situated in one location to ensure that the scale can provide longevity.

Optional Bluetooth available

The M-550 comes with optional Bluetooth capabilities. With the addition of a simple module and a set of protocols, the M-550 is able to integrate into your system with ease.

With a growing focus on telehealth and healthcare in the community, it is important for scales to have the ability to integrate into EMR systems with ease, to centralise data and seamlessly transfer patient information.

Easy to use

The M-550 is a simple and easy to use weighing scale, which features from 4 buttons and a large, clear and easy to ready display.

The On/Off button ensures that you have full control over the device. Hold key enables you to gain weights even for the patient is not steady on the scale and the tare key allows you to deduct any unwanted weight.

Capacity:	160kg	Power Supply:	4 x AA Batteries
Graduations:	200g<100kg>500g	Battery Life:	Over 100 hours of continuous use
Functions:	Hold, Portable, Tare	Approvals:	Class III and MDD Approved
Dimensions:	355mm x 345mm x 60mm	Optional Extras:	Adaptor
Weight of Scale:	4kg	Warranty Period:	4 Years