

MARSDEN

M-560 Premium Slimming and Diet Scale



- 39st capacity
- Accurate to 0.5lb
- Heavy duty construction
- Ideal for outpatient departments
- Suitable for regular transportation
- Optional carry case available

Marsden M-560:

Premium Slimming and Diet Scale

- 39st capacity
- Accurate to 0.5lb
- Perfect for slimming clubs and weight management groups
- Powered by rechargeable battery
- 55 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Optional carry case available to aid portability
- Easy to read remote display
- Hold and Tare functions to enable accurate and easy weighing
- Made from reinforced aluminium structure
- Covered by the Marsden 8 year warranty



Simple weighing

The M-560 is a premium stand on scale, made from reinforced aluminium, which is ideal for weighing adults as it is accurate to a minimum of 0.5lb.

The scale has proven to be a best-seller at Marsden over the years since we first brought this product to market and that is due to the simplicity of the scale and the easy to use functions.

Remote display as standard

The M-560 comes with a remote display fitted as standard. This means that when you are operating the scale in sensitive environments, the weight reading is discreet for only the user to see.

The remote display also means that the user does not have to bend down and see the weight reading on the display, as it is easy to see on the device.

Easy to use

The M-560 uses our popular DP-3810 indicator, which features a full numeric keypad, easy to use functions and a clear display.

The On/Off button ensures that you have full control over the device. Hold enables you to gain weights even for the most active babies and the tare keys allows you deduct the weight of the sling or any clothing.

Capacity:	39st	Power Supply:	Mains/Rechargeable
Graduations:	0.5lb	Battery Life:	55 hours of continuous use
Functions:	Hold, Portable, Tare	Approvals:	Non-Approved
Dimensions:	360mm x 310mm x 90mm	Optional Extras:	Bluetooth/Wi-Fi, Printer, Carry Case
Weight of Scale:	6kg	Warranty Period:	8 Years