

MARSDEN

M-605 Large Hoist Weighing Attachment



- Class III and MDD Approved
- Accurate to 100g
- Larger than the M-600 to suit larger hoists
- Attaches to coat hanger style spreader bars
- Full numeric keypad to calculate BMI
- Optional Wi-Fi or Bluetooth available

Marsden M-605:

Large Hoist Weighing Attachment

- Class III and MDD Approved
- Accurate to 100g
- Attaches to coat hanger style weighing attachments
- Lightweight and easily portable
- Optional carry case available
- Powered by rechargeable battery pack
- 40 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Hold and Tare functions to enable accurate and easy weighing
- Optional Bluetooth and Wi-Fi available
- Covered by the Marsden 8 year warranty



Ideal for care homes

The Marsden M-605 is a hoist weighing attachment which is perfect to weigh immobile patients in care homes, hospitals and private residences.

Simple hook the M-605 to your existing hoist with a coat-hanger style attachments and you are ready to begin gaining accurate readings for your patients and residents. Simple and easy to use.

Optional connectivity options

The M-605 is available with either Bluetooth or Wi-Fi connectivity to ensure that our customers are able to transfer data seamlessly or input into EPR systems.

We would simply need to insert a connectivity module at the time of order and then provide you with our standard protocols for easy and simple integration into your existing system.

Easy to use

The M-605 features from our best-selling DP-3810 indicator. This device features from a full numeric keypad for easy BMI calculation and simple function keys.

When attaching the scale to the hoist, ensure the links are vertical and within a $\pm 5^\circ$ angle to obtain an accurate weight reading. For more information contact the Marsden team.

Capacity:	200kg	Power Supply:	Mains/Rechargeable
Graduations:	100g	Battery Life:	40 hours battery life
Functions:	Hold, Portable, Tare, BMI	Approvals:	Class III and MDD Approved
Dimensions:	730mm x 175mm x 130mm	Optional Extras:	Bluetooth/Wi-Fi, Printer, Carry Case
Weight of Scale:	6kg	Warranty Period:	8 Years