

MARSDEN

M-610 Portable Wheelchair Weigh Beams



- Class III and MDD Approved
- Accurate to 100g
- Tough, robust and hard wearing base
- Suitable for health visitors and hospitals
- Optional carry case available
- Optional Wi-Fi or Bluetooth available

Marsden M-610:

Portable Wheelchair Weigh Beams

- Class III and MDD Approved
- Accurate to 100g
- Perfect for use in the community and hospital
- Powered by rechargeable battery pack
- 40 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Easy to read display
- Hold and Tare functions to enable accurate and easy weighing
- Optional carry case available
- Optional Bluetooth and Wi-Fi available
- Covered by the Marsden 8 year warranty



Simple weighing

The M-610 is a simple and easy to use pair of wheelchair weigh beams, constructed to the highest quality to ensure that it is suitable for weighing patients up to 300kg.

The scale has been designed to be hard wearing making it ideal for heavy use in busy traffic areas and for use out in the community by health visitors who regularly transport their scales.

Optional connectivity options

The M-610 is available with either Bluetooth or Wi-Fi connectivity to ensure that our customers are able to transfer data seamlessly or input into EPR systems.

We would simply need to insert a connectivity module at the time of order and then provide you with our standard protocols for easy and simple integration into your existing system.

Easy to use

The M-610 features from our best-selling DP-3810 indicator. This device features from a full numeric keypad for easy BMI calculation and simple function keys.

The On/Off button ensures that you have full control over the device and you can prolong the battery life. Hold enables you to gain readings for unsteady patients and Tare allows you to deduct unwanted weight.

Capacity:	300kg	Power Supply:	Mains/Rechargeable
Graduations:	100g	Battery Life:	40 hours of continuous use
Functions:	Hold, Portable, Tare	Approvals:	Class III Approved
Dimensions:	960mm x 200mm x 100mm	Optional Extras:	Bluetooth/Wi-Fi, Carry Case, Printer
Weight of Scale:	8.5kg	Warranty Period:	8 Years