

MARSDEN

M-620 Folding Wheelchair Platform Scale



- Class III Approved
- Accurate to 100g
- Wheelchair platform, which can fold for storage
- Wheels as standard for manoeuvrability
- Three displays to show weight, height and BMI
- Optional Wi-Fi or Bluetooth available

Marsden M-620:

Portable, Folding Wheelchair Platform Scale

- Class III Approved
- Accurate to 100g
- Folding base to ensure easy storage
- Wheels fitted as standard for manoeuvrability
- Powered by rechargeable battery pack
- 40 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Easy to read display
- Hold and Tare functions to enable accurate and easy weighing
- Optional Bluetooth and Wi-Fi available
- Covered by the Marsden 8 year warranty



Folding design

The Marsden M-620 has been designed to allow the weighing of almost any wheelchair, motorised or manual, whilst having the ability to be folded and stored efficiently.

With the inclusion of two wheels, the M-620 can be easily folded and wheeled either to a communal storage area and then easily manoeuvred from room to room. Ideal for hospitals, GP surgeries and care homes.

Optional connectivity options

The M-620 is available with either Bluetooth or Wi-Fi connectivity to ensure that our customers are able to transfer data seamlessly or input into EPR systems.

We would simply need to insert a connectivity module at the time of order and then provide you with our standard protocols for easy and simple integration into your existing system.

Easy to use

The M-620 features from our best-selling DP-3810 indicator. This device features from a full numeric keypad for easy BMI calculation and simple function keys.

The On/Off button ensures that you have full control over the device and you can prolong the battery life. Hold enables you to gain readings for unsteady patients and Tare allows you to deduct unwanted weight.

Capacity:	300kg	Power Supply:	Mains/Rechargeable
Graduations:	100g	Battery Life:	40 hours of continuous use
Functions:	Hold, Portable, Tare, BMI	Approvals:	Class III Approved
Dimensions:	740mm x 700mm x 1160mm	Optional Extras:	Bluetooth/Wi-Fi, Printer
Weight of Scale:	38kg	Warranty Period:	8 Years