

# MARSDEN

## M-640 Professional Wheelchair Scale



- Class III Approved
- Accurate to 50g
- Multi-functional wheelchair weigher
- Seat for patients unstable on their feet
- Ramps for wheelchair bound patients
- Optional Wi-Fi or Bluetooth available

# Marsden M-640:

## Professional, Multi-Functional Wheelchair Scale

- Class III Approved
- Accurate to 50g
- Multi-functional wheelchair scale
- Folding seat for patients unsteady on their feet
- Gradual incline ramps for wheelchair bound patients
- Powered by rechargeable battery pack
- 55 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Hold and Tare functions to enable accurate and easy weighing
- Optional Bluetooth and Wi-Fi available
- Covered by the Marsden 4 year warranty



### Multi-functional use

The M-640 is a high capacity medical wheelchair scale with a convenient design to weigh bariatric patients and those who are unable to stand unassisted.

The extra large base is able to accommodate almost any wheelchair scale. The chair is ideal for elderly and patients which are unsteady on their feet and sturdy handrails are suitable for all.

### Optional connectivity options

The M-640 is available with either Bluetooth or Wi-Fi connectivity to ensure that our customers are able to transfer data seamlessly or input into EPR systems.

We would simply need to insert a connectivity module at the time of order and then provide you with our standard protocols for easy and simple integration into your existing system.

### Easy to use

The M-640 features from our best-selling DP-3810 indicator. This device features from a full numeric keypad for easy BMI calculation and simple function keys.

The On/Off button ensures that you have full control over the device and you can prolong the battery life. Hold enables you to gain readings for unsteady patients and Tare allows you to deduct unwanted weight.

|                         |                                   |                         |                          |
|-------------------------|-----------------------------------|-------------------------|--------------------------|
| <b>Capacity:</b>        | 300kg or 500g                     | <b>Power Supply:</b>    | Mains/Rechargeable       |
| <b>Graduations:</b>     | 50g<150kg>100g or 100g<300kg>200g | <b>Battery Life:</b>    | 55 hours of battery life |
| <b>Functions:</b>       | Hold, Portable, Tare, BMI         | <b>Approvals:</b>       | Class III Approved       |
| <b>Dimensions:</b>      | 800mm x 800mm                     | <b>Optional Extras:</b> | Bluetooth/Wi-Fi, Printer |
| <b>Weight of Scale:</b> | 53kg                              | <b>Warranty Period:</b> | 4 Years                  |