

MARSDEN

M-653 Wheelchair Scale with Folding Handrail



- Class III and MDD Approved
- Accurate to 100g
- Portable, wheelchair weigher
- Handrail to give support to unsteady patients
- Gentle incline ramps
- Optional Wi-Fi or Bluetooth available

Marsden M-653:

Portable Wheelchair Scale with Folding Handrail

- Class III and MDD Approved
- Accurate to 100g
- Portable, wheelchair weigher
- Handrail to help support patients unsteady on their feet
- Gradual incline ramps
- Powered by rechargeable battery pack
- 55 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Hold and Tare functions to enable accurate and easy weighing
- Optional Bluetooth and Wi-Fi available
- Covered by the Marsden 4 year warranty



Handrail to support

The Marsden M-653 is our standard best-selling M-650 platform, with an additional handrail to provide support for patients which are unsteady on their feet.

The handrail on the M-653 also has a hinge, which means that the rail can easily be folded down and then safely manoeuvred into position, or conveniently stored easily thanks to the handle and wheels.

Optional connectivity options

The M-653 is available with either Bluetooth or Wi-Fi connectivity to ensure that our customers are able to transfer data seamlessly or input into EPR systems.

We would simply need to insert a connectivity module at the time of order and then provide you with our standard protocols for easy and simple integration into your existing system.

Easy to use

The M-65 features from our best-selling DP-3810 indicator. This device features from a full numeric keypad for easy BMI calculation and simple function keys.

The On/Off button ensures that you have full control over the device and you can prolong the battery life. Hold enables you to gain readings for unsteady patients and Tare allows you to deduct unwanted weight.

Capacity:	300kg	Power Supply:	Mains/Rechargeable
Graduations:	100g	Battery Life:	55 hours of continuous use
Functions:	Hold, Portable, Tare, BMI	Approvals:	Class III Approved
Dimensions:	1150mm x 800mm x 66mm	Optional Extras:	Bluetooth/Wi-Fi, Printer
Weight of Scale:	28.6kg	Warranty Period:	4 Years