

# MARSDEN

## M-800 Bariatric and Geriatric Scale



- Class III Approved
- Accurate to 50g
- Large platform with built-in handrails
- Heavy duty weighing scale
- Suitable for bariatric patients
- Optional Wi-Fi or Bluetooth available

# Marsden M-800:

## Bariatric and Geriatric Weighing Scale

- Class III Approved
- Accurate to 50g
- Large platform for patients to stand on
- Sturdy handrails to provide support
- Heavy duty construction
- Powered by rechargeable battery pack
- 55 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Hold and Tare functions to enable accurate and easy weighing
- Optional Bluetooth and Wi-Fi available
- Covered by the Marsden 4 year warranty



### Robust weighing scale

The Marsden M-800 is a high capacity weighing scales with a large weighing platform and sturdy handrails to give patients support.

This scale has been designed with geriatric and bariatric patients in mind. The low profile platform means that patients can easily step on the scale, making the process easier and safer.

### Optional connectivity options

The M-800 is available with either Bluetooth or Wi-Fi connectivity to ensure that our customers are able to transfer data seamlessly or input into EPR systems.

We would simply need to insert a connectivity module at the time of order and then provide you with our standard protocols for easy and simple integration into your existing system.

### Easy to use

The M-800 features from our best-selling DP-3810 indicator. This device features from a full numeric keypad for easy BMI calculation and simple function keys.

The On/Off button ensures that you have full control over the device and you can prolong the battery life. Hold enables you to gain readings for unsteady patients and Tare allows you to deduct unwanted weight.

<b>Capacity:</b>	300kg	<b>Power Supply:</b>	Mains/Rechargeable
<b>Graduations:</b>	50g<150kg>100g	<b>Battery Life:</b>	55 hours of continuous use
<b>Functions:</b>	Hold, Portable, Tare, BMI	<b>Approvals:</b>	Class III Approved
<b>Dimensions:</b>	700mm x 600mm x 50mm	<b>Optional Extras:</b>	Bluetooth/Wi-Fi, Printer
<b>Weight of Scale:</b>	33kg	<b>Warranty Period:</b>	4 Years