

MARSDEN

M-900 Portable Bed Weighing Scales



- Class III Approved
- Accurate to 200g
- Tough and robust
- Suitable for weighing in hospital and care homes
- Low profile ramps
- Optional Wi-Fi or Bluetooth available

Marsden M-900:

Portable Bed Weighing Beams

- Class III Approved
- Accurate to 200g
- Perfect for use in hospitals and care homes
- Powered by rechargeable battery pack
- 55 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Low profile ramps to ensure ease when pushing a bed
- Easy to read display
- Hold and Tare functions to enable accurate and easy weighing
- Optional Bluetooth and Wi-Fi available
- Covered by the Marsden 4 year warranty



Designed for ease of use

The M-900 was designed to help make weighing immobile patients in their bed a much easier and stress-free process for healthcare professionals.

Place the weigh beams behind the castors of the bed and then slowly and carefully push the bed up the ramps. The weight will then show on the display and Tare can be used to deduct the weight of the bed.

Optional connectivity options

The M-900 is available with either Bluetooth or Wi-Fi connectivity to ensure that our customers are able to transfer data seamlessly or input into EPR systems.

We would simply need to insert a connectivity module at the time of order and then provide you with our standard protocols for easy and simple integration into your existing system.

Easy to use

The M-900 features from our best-selling DP-3810 indicator. This device features from a full numeric keypad for easy BMI calculation and simple function keys.

The On/Off button ensures that you have full control over the device and you can prolong the battery life. Hold enables you to gain readings for unsteady patients and Tare allows you to deduct unwanted weight.

Capacity:	600kg	Power Supply:	Mains/Rechargeable
Graduations:	200g	Battery Life:	55 hours of continuous use
Functions:	Hold, Portable, Tare, BMI	Approvals:	Class III Approved
Dimensions:	1240mm x 360mm x 65mm	Optional Extras:	Bluetooth/Wi-Fi, Printer
Weight of Scale:	11.5kg	Warranty Period:	4 Years