


MARSDEN

M-910 Floor Mounted Bed Scale with Ramps




- Class III Approved
- Accurate from 200g
- Tough and robust weighing platform
- Suitable for weighing beds in hospitals
- Low profile ramps for patient comfort
- Optional Wi-Fi or Bluetooth available

Marsden M-910:

Bed and Trolley Weigher

- Class III Approved
- Accurate from 200g
- Perfect for weighing beds and trolleys in hospitals
- Powered by rechargeable battery pack
- 55 hours of continuous use
- Auto power-off function helps to save battery power
- Simple and easy to use
- Low profile ramps to ensure ease when pushing a bed
- Easy to read display
- Hold and Tare functions to enable accurate and easy weighing
- Optional Bluetooth and Wi-Fi available
- Covered by the Marsden 4 year warranty



Designed for ease of use

The M-910 is a bed weighing platform scale. Built to ensure that it is easy for medical professionals to safely and accurately obtain a weight for their patients.

This scale features from Tare function which enables staff to remove the weight of the bed or trolley, to ensure that they are left with an accurate weight for their patient, without having to calculate the difference.

Optional connectivity options

The M-910 is available with either Bluetooth or Wi-Fi connectivity to ensure that our customers are able to transfer data seamlessly or input into EPR systems.

We would simply need to insert a connectivity module at the time of order and then provide you with our standard protocols for easy and simple integration into your existing system.

Easy to use

The M-910 features from our best-selling DP-3810 indicator. This device features from a full numeric keypad for easy BMI calculation and simple function keys.

The On/Off button ensures that you have full control over the device and you can prolong the battery life. Hold enables you to gain readings for unsteady patients and Tare allows you to deduct unwanted weight.

Capacity:	600kg or 1500kg	Power Supply:	Mains/Rechargeable
Graduations:	200g or 500g	Battery Life:	55 hours of continuous use
Functions:	Hold, Tare, BMI	Approvals:	Class III Approved
Dimensions:	1200mm x 1000mm / 1200mm x 1200mm / 1500mm x 1500mm	Warranty Period:	4 Years