

# MARSDEN

## M-950 Portable Bed Weigher



- Class III Approved
- Accurate from 200g
- Four pad weighing system
- Suitable for beds in hospitals and care homes
- Can be used with most patient beds
- Optional Wi-Fi or Bluetooth available

# Marsden M-950:

## Portable Bed Weighing System

- Class III Approved
- Accurate from 200g
- Perfect for weighing beds and trolleys in hospitals and care homes
- Powered by rechargeable battery pack
- 40 hours of continuous use
- Auto power-off function helps to save battery power
- Four pad weighing system which can be stored on our trolley
- Shallow ramps to ensure patient comfort
- Easy to read display
- Hold and Tare functions to enable accurate and easy weighing
- Optional Bluetooth and Wi-Fi available
- Covered by the Marsden 8 year warranty



### Why the M-950?

The M-950 is an easy to use, portable bed weighing scale. It makes weighing bed-bound patients straightforward and completely hassle free with no discomfort.

This scale consists of four weighing pads. Simply place the pads behind the castors, gently push the bed up the shallow ramps and the weight will appear on the screen. Tare can deduct the weight of the bed.

### Optional connectivity options

The M-950 is available with either Bluetooth or Wi-Fi connectivity to ensure that our customers are able to transfer data seamlessly or input into EPR systems.

We would simply need to insert a connectivity module at the time of order and then provide you with our standard protocols for easy and simple integration into your existing system.

### Easy to use

The M-950 features from our best-selling DP-3810 indicator. This device features from a full numeric keypad for easy BMI calculation and simple function keys.

The On/Off button ensures that you have full control over the device and you can prolong the battery life. Hold enables you to gain readings for unsteady patients and Tare allows you to deduct unwanted weight.

<b>Capacity:</b>	600kg or 1000kg	<b>Power Supply:</b>	Mains/Rechargeable
<b>Graduations:</b>	200g or 500g	<b>Battery Life:</b>	40 hours of continuous use
<b>Functions:</b>	Hold, Tare, BMI, Portable	<b>Approvals:</b>	Class III Approved
<b>Dimensions:</b>	650mm x 480mm x 980mm	<b>Optional Extras:</b>	Bluetooth/Wi-Fi
<b>Weight Per Pad:</b>	5kg each	<b>Warranty Period:</b>	8 Years