

MARSDEN

M-955 Large Bed Weigher



- Class III Approved
- Accurate from 200g
- Four pad weighing system for large beds
- Suitable for beds in hospitals and care homes
- Can be used with most patient beds
- Optional Wi-Fi or Bluetooth available

Marsden M-950:

Large Bed Weighing System

- Class III Approved
- Accurate from 200g
- Perfect for weighing beds and trolleys in hospitals and care homes
- Powered by rechargeable battery pack
- 40 hours of continuous use
- Auto power-off function helps to save battery power
- Four pad weighing system which can be stored on our trolley
- Designed for use with beds featuring double castors
- Easy to read display
- Hold and Tare functions to enable accurate and easy weighing
- Optional Bluetooth and Wi-Fi available
- Covered by the Marsden 8 year warranty



Why the M-955?

The M-955 is a high capacity, Class III Approved bed weighing scale, designed to weigh large hospital beds which feature from double castors.

These weigh pads are slightly longer than our standard M-950, which enables the M-955 to accommodate a wider variety of beds, including Medstrom brand. With the same benefits as the standard M-950.

Optional connectivity options

The M-955 is available with either Bluetooth or Wi-Fi connectivity to ensure that our customers are able to transfer data seamlessly or input into EPR systems.

We would simply need to insert a connectivity module at the time of order and then provide you with our standard protocols for easy and simple integration into your existing system.

Easy to use

The M-955 features from our best-selling DP-3810 indicator. This device features from a full numeric keypad for easy BMI calculation and simple function keys.

The On/Off button ensures that you have full control over the device and you can prolong the battery life. Hold enables you to gain readings for unsteady patients and Tare allows you to deduct unwanted weight.

Capacity:	600kg or 1000kg	Power Supply:	Mains/Rechargeable
Graduations:	200g or 500g	Battery Life:	40 hours of continuous use
Functions:	Hold, Tare, BMI, Portable	Approvals:	Class III Approved
Dimensions:	650mm x 480mm x 980mm	Optional Extras:	N/A
Weight Per Pad:	10kg each	Warranty Period:	8 Years