

# MARSDEN

## MBF-6000 Body Composition Scale



- Class III and MDD Approved
- Accurate to 100g
- Provides a complete picture of your body health
- Wide variety of readings provided
- Closer to DEXA than any body composition scale
- Sold with a printer as standard

# Marsden MBF-6000:

## Body Composition Scale

- Class III and MDD Approved
- Accurate to 100g
- Sold with a printer as standard
- Three weight displays for weight, height and BMI
- 55 hours of continuous use
- Provides various body composition readings
- Auto power-off function helps to save battery power
- Simple and easy to use
- Closer to DEXA than any other body composition scales
- Easy to read display
- Hold and Tare functions to enable accurate and easy weighing
- Covered by the Marsden 4 year warranty



### The readings

The Marsden MBF-6000 is a body composition weighing scale, which is ideal for use in hospitals, gyms, personal trainers and sports clubs.

Measures Basal Metabolic Rate, Total Body Water), Fat Free Mass, Fat Mass, Body Fat %, Muscle Mass, Protein Mass, Extracellular Water, Intracellular Water, Skeletal Muscle , Visceral Fat Area Level, Metabolic Age (AGEM), Weight and BMI

### Closest to DEXA

This scale was used in a study by Hull/ York Medical School, who compared body composition scales and found that the MBF was closer to DEXA than any other.

DEXA is considered to the gold standard, when it comes to body composition analysers. So much so, that it is used by the medical profession and top athletes, to monitor themselves and ensure that they are in peak performance.

### Easy to use

The MBF-6000 is a simple and easy to use body composition analyser. Thanks to the full numeric indicator and two displays, everything is easy to input and see.

When your patient steps onto this scale, a very safe and small electrical current passes through their body. Electrical current flows through water, fat and muscle mass at different speeds, and therefore the scale can calculate how a body is composed.

<b>Capacity:</b>	300kg	<b>Power Supply:</b>	Mains/Rechargeable
<b>Graduations:</b>	100g	<b>Battery Life:</b>	55 hours of continuous use
<b>Functions:</b>	Hold, Portable, Tare, BMI	<b>Approvals:</b>	Class III and MDD Approved
<b>Dimensions:</b>	450mm x 340mm x 85mm	<b>Optional Extras:</b>	Carry Case
<b>Weight of Scale:</b>	9.5kg	<b>Warranty Period:</b>	4 Years